# Imunoglobulina A

Imunoglobulina A / Imunoglobulina A  
Ref. 20.010.00

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### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

### INTENDED USE

Kit intended to determine Immunoglobulin A (IgA) in serum and plasma samples. Diagnostic use only.

### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

### WORKING PRINCIPLE

**Method:** Immunoturbidimetric.

Anti-immunoglobulin A antibodies form insoluble complexes when mixed with samples containing IgA. The turbidity intensity is proportional to the IgA concentration of the sample.

### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum and plasma (heparin or EDTA).

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum / Plasma	4 to 8 °C	8 months
	-20 °C	8 months

### PRODUCT DESCRIPTION

**R 1**

Goat anti-Immunoglobulin A antibodies; activator; phosphate buffer; preservative.

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and control below: Multicalibrador Plus **REF** 21.011.00 Medium Multicalibrador Plus **REF** 21.003.00

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Test Tubes, glass pipettes and/or automatic.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

**Automation:** This product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

<b>Reaction Type</b>	End Point
<b>Wave length</b>	340 nm
<b>Temperature</b>	37 °C
<b>Sample Volume *</b>	2 µL
<b>R1 Volume *</b>	285 µL
<b>R1 Incubation Time</b>	180 seconds
<b>R1 + Sample Reaction Time</b>	120 seconds
<b>Calibration **</b>	Calibration curve

\* Sample and R1 volumes can be modified as long as the established ratio is maintained.

### \*\* CALIBRATION CURVE:

Dilute the Multicalibrador Plus in a 9,0 g/L NaCl solution as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	14	10	25	50	100
NaCl 9,0 g/L (µL)	100	386	90	75	50	-
Dilution Factor	0	0,035	0,1	0,25	0,5	1

Multiply the Immunoglobulin A concentration indicated in the Multicalibrador Plus instruction of use by the corresponding dilution factor to obtain the IgA concentration of each dilution. The IgA concentration of the Multicalibrador Plus is traceable to the Standard Reference Material ERM DA470k / IFCC.

### B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

### C) INTERPRETATION

Quantification of serum immunoglobulins is important for the diagnosis of primary or secondary immunodeficiencies, for the monitoring of immunoglobulin therapy, and for the follow-up of the clinical course of multiple myeloma. IgA accounts for about 10 to 15% of the total circulating immunoglobulins. The most important form of IgA is secretory, found in saliva, tears, sweat, colostrum, milk, bronchial and gastrointestinal secretions. The decrease in IgA levels occur mainly due to congenital and acquired immunodeficiencies. Elevation of polyclonal IgA is a response to infections of the skin, lungs, kidneys and cirrhosis of the liver. Increases in monoclonal IgA (paraprotein) are found in multiple myeloma and in plasma cell proliferative disorders.

### PERFORMANCE CHARACTERISTICS

Operating range
19 to 601 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1.17 mg/dL	18.73 mg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
100 mg/dL	20 mg/dL	1500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 1.0323x – 12.515
Correlation Coefficient (R)	0.9978

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 80 mg/dL, 180 mg/dL and 500 mg/dL was -12.4%, -3.7% and 0.7%, respectively.

### Prozone High Dose Effect:

Prozone effect does not occur up to 1060.88 mg/dL of Immunoglobulin A.

### Precision:

The intra-assay precision studies were performed in 5 replicates for 10 days; the inter-assay studies were performed in 3 conditions, with 3 replicates for 5 days.

Samples (mg/dL)	Intra-assay Precision		Samples (mg/dL)	Inter-assay Precision	
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)
95.54	1.53	1.6	96.70	1.32	1.4
187.87	2.95	1.6	191.88	4.35	2.3
405.80	9.16	2.3	409.55	10.31	2.5

CV: Coefficient of variation; SD: Standard deviation

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

### REFERENCE RANGES

Newborn (4 days)	0 a 2.2 mg/dL
Adults (20-60 years)	70 a 400 mg/dL
Adults (> 60 anos)	90 a 410 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling for +55 35 3214 4646.

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).

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- PRICE, C. P.; et al. Ann Clin Biochem. v. 20, p. 1-14, 1983.
- DATI, F. et al; Eur J Clin Chem Clin Biochem. v. 34, p. 517-520, 1996.
- TIETZ, Textbook of Clinical Chemistry. 3 ed. Burtis CA, Ashwood ER. WB Saunders Co, 1999.
- FRIEDMAN AND YOUNG, Effects of the disease on clinical laboratory tests, 3 ed., AACCC Press, 1997.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACCC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.BURTIS, C.A.; ASHWOOD, E.R.; BRUNS, D.E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.

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