Imunoglobulina M

Imunoglobulina M / Imunoglobulina M
Ref. 20.012.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine immunoglobulin M (IgM) in serum and plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Immunoturbidimetric.

Anti-human immunoglobulin M antibodies form insoluble complexes when mixed with samples containing immunoglobulin M. The turbidity intensity is proportional to the immunoglobulin M concentration in the sample, which is spectrophotometrically measured at 340 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (heparin or EDTA).

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

| | Temperature | Sample Stability |
|------------------|-------------|------------------|
| Serum and Plasma | 4 to 8 °C | 4 months |
| | -20 °C | 6 months |

PRODUCT DESCRIPTION

R 1

Goat anti-human immunoglobulin M antibodies; activator; phosphate buffer; preservative.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and control below:

| | | |
|-----------------------------|------------|-----------|
| Multicalibrador Plus | REF | 21.011.00 |
| Medium Multicalibrador Plus | | 21.003.00 |

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Test Tubes, glass pipettes and/or automatic.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Automation: This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaaltda.ind.br

| Reaction Type | End Point |
|---------------------------|-------------------|
| Wavelength | 340 nm |
| Temperature | 37 °C |
| Sample Volume / CAL * | 3,0 µL |
| R1 Volume * | 300 µL |
| R1 Incubation Time | 180 seconds |
| R1 + Sample Reaction Time | 180 seconds |
| Calibration ** | Calibration curve |

* Sample and R1 volumes can be modified as long as the established ratio is maintained.

** CALIBRATION CURVE:

Dilute the Multicalibrador Plus in a 9.0 g/L NaCl solution as follows:

| Dilution | 1 | 2 | 3 | 4 | 5 | 6 |
|-------------------|-----|------|------|------|------|-----|
| Calibrador (µL) | - | 10 | 25 | 50 | 75 | 100 |
| NaCl 9.0 g/L (µL) | 100 | 90 | 75 | 50 | 25 | - |
| Dilution Factor | 0 | 0.10 | 0.25 | 0.50 | 0.75 | 1 |

Multiply the immunoglobulin M concentration indicated in the Multicalibrador Plus instruction of use by the corresponding dilution factor to obtain the immunoglobulin M concentration of each dilution. The immunoglobulin M concentration of the Multicalibrador Plus is traceable to the Standard Reference Material ERM DA470K/IFCC.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

Quantification of serum immunoglobulins is important for the diagnosis of primary or secondary immunodeficiencies, for the monitoring of immunoglobulin therapy, and for the follow-up of the clinical course of multiple myeloma. Immunoglobulin M is the only immunoglobulin synthesized by newborns. It accounts for about 5% to 10% of total circulating immunoglobulins. It is a pentameric molecule, and its large size prevents its passage into extravascular spaces. Its decrease occurs in acquired or congenital immunodeficiencies. Elevation of polyclonal immunoglobulin M (normal response to infections) occurs mainly in primary viral infections and infections of the bloodstream, such as malaria. It may also be increased in primary biliary cirrhosis and chronic active hepatitis. Increases in monoclonal immunoglobulin M (paraprotein) are found in plasma cell proliferative disorders, such as Waldenström macroglobulinemia.

PERFORMANCE CHARACTERISTICS

| Operating range |
|----------------------|
| 15,36 a 303,04 mg/dL |

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

| Sensitivity | |
|-----------------|----------------------|
| Detection Limit | Quantification Limit |
| 0,83 mg/dL | 15,36 mg/dL |

| Analytical Specificity | | | |
|------------------------|-----------|---------------|-------------------|
| Hemoglobin | Bilirubin | Triglycerides | Rheumatoid factor |
| 100 mg/dL | 20 mg/dL | 2000 mg/dL | 150 UI/mL |

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy | |
|-----------------------------|----------------------|
| Number of Samples | 35 in duplicate |
| Regression Equation | y = 0.9464x + 0.2748 |
| Correlation Coefficient (R) | 0.9934 |

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 40 mg/dL, 230 mg/dL and 300 mg/dL was -4.70%, -5.20% and -5.30%, respectively.

Prozone High Dose Effect:

Prozone effect does not occur up to 549.01 mg/dL of immunoglobulin M.

Precision:

The intra-assay precision studies were performed in 5 replicates for 10 days; the inter-assay studies were performed in 3 conditions, with 3 replicates for 5 days.

| Samples (mg/dL) | Intra-assay Precision | | Samples (mg/dL) | Inter-assay Precision | |
|-----------------|-----------------------|--------|-----------------|-----------------------|--------|
| | SD (mg/dL) | CV (%) | | SD (mg/dL) | CV (%) |
| 63.75 | 0.79 | 1.20 | 62.55 | 0.80 | 1.30 |
| 146.27 | 1.09 | 0.70 | 145.70 | 1.03 | 0.70 |
| 261.77 | 2.34 | 0.90 | 260.88 | 2.41 | 0.90 |

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

| | |
|------------------------------|----------------|
| Newborns (4 days) | 5 - 30 mg/dL |
| Adults (20 to 60 years) | 40 - 230 mg/dL |
| Adults (older than 60 years) | 30 - 360 mg/dL |

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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