Canines	500 mg/dL	40 mg/dL	1000 mg/dL
Felines	500 mg/dL	40 mg/dL	1000 mg/dL
Equines	500 mg/dL	5 mg/dL	500 mg/dL
Bovines	500 mg/dL	5 mg/dL	500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy – Canines	
Number of Samples	40 in duplicate
Regression Equation	y = 1.0082x – 0.3798
Correlation Coefficient (R)	0.9988
Accuracy – Felines	
Number of Samples	34 in duplicate
Regression Equation	y = 0.9536x + 6.2541
Correlation Coefficient (R)	0.9909
Accuracy – Equines	
Number of Samples	16 in duplicate
Regression Equation	y = 0.9868x + 1.9084
Correlation Coefficient (R)	0.9990
Accuracy – Bovines	
Number of Samples	16 in duplicate
Regression Equation	y = 1.0184x – 0.7314
Correlation Coefficient (R)	0.9967

#### Precision:

Determined with two runs.

Samples (mg/dL)	Within-Run Precision		Samples (mg/dL)	Total Precision	
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)
116.95	2.71	2.3	117.53	2.41	2.1
317.64	10.15	3.2	310.69	10.35	3.3

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Canines	135 – 270 mg/dL
Felines	95 – 130 mg/dL
Equines	75 – 150 mg/dL
Bovines	80 – 120 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units (mmol/L):

Cholesterol (mg/dL) x 0,026 = Cholesterol (mmol/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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#### Conversion to the Unit of the International System (mmol/L):

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- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. **Clinical biochemistry of domestic animals**. 5th. ed. San Diego: Academic Press, c1997 932p.
- BELL, L. L. et al. A simplified method for the estimation of total cholesterol in serum and demonstration of its specificity. **J. Biol. Chem.** v.195, p.357-366, 1952.
- TRINDER, P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. **Ann. Clin. Biochem.** v.6, p.24-27, 1969.
- BURTIS, C. A.; ASHWOOD, E. R.; BRUNS, D. E. **Tietz Fundamentos de Química Clínica**, Saunders Elsevier, 6 ed, 2008.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.

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