	
LDL Colesterol	
LDL Cholesterol / LDL Colesterol Ref. 10.015.00	
Responsável Técnico: Dr. Gilson Serio Pizzo CRF MG – 5310 Anvisa 80027310233	

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine LDL cholesterol in serum and plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use. Lyophilized calibrator.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

PREPARATION AND HANDLING

Calibrator: reconstitute with 1 mL of purified water. Homogenize gently by inversion, avoiding foam formation. Leave to rest for 30 minutes at room temperature until the product is completely dissolved.
 Reconstituted, the calibrator is stable for 14 days when stored at 2 to 8 °C, and 3 months when stored at -20 °C.

WORKING PRINCIPLE

Method: Selective surfactant / Direct
 R1 contains polyvinyl sulphonic acid (PVS) and polyethylene glycol methyl ether (PEGME) which forms a complex with LDL-VLDL-chylomicrons, leaving HDL cholesterol free. This is consumed by Cholesterol oxidase and Cholesterol esterase. After the addition of R2 the LDL is released from the PVS / PEGME complex by specific detergents. Free LDL reacts with the enzymes producing hydrogen peroxide which is quantified by the Trinder reaction. The intensity of the color produced is proportional to the amount of LDL cholesterol present in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum or plasma (heparin and EDTA)
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

	Temperature	Stability Period
Serum and Plasm	4 to 8 °C	7 days
	-20 °C	2 months

PRODUCT DESCRIPTION

R 1	MES buffer pH 6.5; polyvinylsulfonic acid, polyethylene glycol methyl ester, Magnesium Chloride, Detergents, EDTA, 4-aminoantipyrine, Cholesterol esterase, Cholesterol Oxidase and Peroxidase.
	MES buffer pH 6.5, EDTA, TODB N Detergent, N-bis (4-sulfobutyl) -3-methylaniline
R 2	Human serum, preservatives. Value printed on the bottle label.
	Traceable to CDC reference method for LDL Cholesterol.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For the Internal Quality Control Laboratory it is recommended:

LDL Cholesterol CAL (included in kit)	10.015.00
Normal Control - Quantinorm	13.003.00
Pathological Control - Quantialt	13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 600 nm.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.tda.ind.br

Types of reaction	Final/Growing Point
Wave length	600 nm
Temperature	37 °C
Sample volume, blank (Purified water) or calibrator *	3 µL
Volume R1*	225 µL
Reading 1 (A1)	5 minutes
Volume R2*	75 µL
Reading 2 (A2)	5 minutes
Calibration	Linear

* Volumes can be modified as long as the established ratio is maintained.

B) CALCULATIONS

$$\text{LDL Cholesterol (mg/dL)} = \frac{\Delta A \text{ Sample} - \Delta A \text{ Blank}}{\Delta A \text{ Calibrador} - \Delta A \text{ Blank}} \times \text{Calibrator Value (mg/dL)}$$

Where:

$$\Delta A = A2 - A1$$

C) INTERPRETATION

Very low density lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on VLDL (Very Low Density Lipoproteins). The elimination of LDL cholesterol from plasma occurs primarily by the liver in parenchymal cells through specific LDL receptors. Most of the cholesterol present in atherosclerotic plaques originates from LDL. For this reason LDL cholesterol is considered an important clinical predictor of coronary atherosclerosis.

PERFORMANCE CHARACTERISTICS,

Operating range
5.50 to 241 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1.50 mg/dL	5.47 mg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
1000 mg/dL	40 mg/dL	1500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum
Number of Samples	37 in duplicate
Regression Equation	y = 0.994x + 2.02
Correlation Coefficient (R)	0.9953

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 52 mg/dL and 220 mg/dL was 3.28 % and 0.32%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (µmol/L)	Replicates	Within-Run Precision		Total Precision	
		SD (µmol/L)	CV (%)	SD (µmol/L)	CV (%)
32,401	80	0,580	1,8	1,019	3,1
62,944	80	0,689	1,1	1,023	1,6
160,263	80	0,947	0,6	1,423	0,9

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Suggested therapeutic target values for cardiovascular risk category estimated by the Doctor for adults over 20 years:

Risk	Therapeutic target (mg/dL)
Low	< 130
Intermediary	< 100
High	< 70
Very High	< 50

NOTE: According to the Brazilian Consensus for the Normalization of Laboratory Determination of the Lipid Profile, reference values for LDL cholesterol were replaced by the values of therapeutic target, stratified according to the calculated cardiovascular risk determination of the Lipid Profile.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- BADIMON JJ, BADIMON L, FUESTER V. Regression of atherosclerotic lesions by high density lipoprotein plasma fraction in the cholesterol-fed rabbit. *Journal of Clinical Investigation*, 85; 1234-41 (1990).
- BARR DP, RUSS EM, EDER HA. Protein-lipid relationships in human plasma. *Am. J. Med.*, 11;480-493 (1951).
- CASTELLI WP et al. HDL Cholesterol and other lipids in coronary heart disease. *Circulation* ahajournals.org. 55;767-772 (1977).
- CROUSE JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease. *J. Lipid Res.* 26, 566-574 (1985).
- FALUDI AA, IZAR MCO, SARAIVA JFK, CHACRA APM, BIANCO HT, AFIUNE NETO A et al. Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose -2017. *Arq Bras Cardiol* 2017; 109(2Supl.1):1-76.
- GORDON T et al., High density lipoprotein as a protective factor against coronary heart disease. *Am. J. Med.*, 62; 707-714 (1977).
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACCPress, 2000.

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