	
LDL Direto Plus	
LDL Direct Plus / LDL Directo Plus Ref. 10.043.00	
Responsável Técnico: Dr. Gilson Serio Pizzo CRF MG – 5310 Anvisa 80027310310	

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE
 Kit intended to determine LDL in serum and plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product in use is stable for 4 weeks when stored *on board*.
- Do not use reagents whose shelf life has expired.

PREPARATION AND HANDLING

Calibrator: reconstitute with 1 mL of purified water. Homogenize gently by inversion, avoiding foam formation. Wait 30 minutes at room temperature until the product is completely dissolved.

Reconstituted, the calibrator is stable for 30 days when stored at 2 to 8 °C or -20 °C.

WORKING PRINCIPLE

Method: Direct
 The polymer and the surfactant present in R1 form a compound with LDL, selectively inhibiting the reaction of LDL with cholesterol esterase and solubilizing chylomicrons and HDL and VLDL lipoproteins. Solubilized cholesterol is consumed by cholesterol esterase and cholesterol oxidase in a colorless reaction. In the second step of the reaction, the LDL particles are solubilized by a specific surfactant present in R2 and the LDL-cholesterol is released to react with cholesterol esterase and cholesterol oxidase forming a colored compound that is spectrophotometrically measured at 546 nm. The intensity of the color produced is proportional to the amount of LDL cholesterol present in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (EDTA and heparin)

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	7 days
	-20 °C	3 months

PRODUCT DESCRIPTION

R 1	Buffer MOPS pH 7.0; 4-aminoantipyrine 0.6 mmol/L; colesterol esterase 2 KU/L; colesterol oxidase 5 KU/L; Peroxidase 15 KU/L; sodium cholate 100 mmol/L; magnesium chloride 10 mmol/L
	Buffer MOPS pH7; TOOS 4 mmol/L; triton X-100 100 mmol/L; magnesium chloride 10 mmol/L
R 2	Buffer MOPS pH7; TOOS 4 mmol/L; triton X-100 100 mmol/L; magnesium chloride 10 mmol/L

CAL	LDL, sodium chloride 100 mmol/L, BSA 20 g/L, sucrose 10% and preservative 0.03%.	
------------	--	---

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Control of Laboratory Quality, it is recommended to use the calibrator and the controls below:

Calibrator – LDL Direct CAL	REF	10.043.00
Normal control – Quantinorm		13.003.00
Pathological control - Quantialt		13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 546 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.ind.br

Reaction Type	Two Point / Positive
Wavelength	546 nm
Temperature	37 °C
Sample volume*	3 µL
R1 Volume*	300 µL
Sample + R1 Incubation Time	300 seconds
First Reading	Immediately
R2 Volume*	100 µL
R1 + Sample + R2 Incubation Time	300 seconds
Second Reading	Immediately
Calibration	2 points

* Volumes can be modified as long as the established proportion is maintained.

B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample. Factor calculation

C) INTERPRETATION

Cholesterol is a lipid that constitutes cell membranes and is a precursor to steroid hormones and bile acids that are synthesized by the body's cells and absorbed with food. The transport of cholesterol in plasma occurs through lipoproteins that are lipids and apolipoproteins complexes. There are four classes of lipoproteins: high density lipoproteins (HDL), low density lipoproteins (LDL), very low density lipoproteins (VLDL) and chylomicrons.

Low density lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on VLDL. The elimination of LDL cholesterol from plasma occurs mainly by the liver in parenchymal cells through specific LDL receptors. Most of the cholesterol present in atherosclerotic plaques originates from LDL. It is important to test LDL cholesterol in addition to total cholesterol, because even with the normal total cholesterol level, an increase in LDL cholesterol can occur with an increased risk of coronary artery disease.

PERFORMANCE CHARACTERISTICS

Operating range
11,54 a 496,15 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Detection Limit
11,54 mg/dL

Analytical Specificity		
Bilirubin	Hemoglobin	Triglycerides
14,6 mg/dL	500 mg/dL	1610 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

The accuracy test was performed with 3 lots by comparison with a reference material.

Accuracy			
Lot	Reference value (mg/dL)	Obtained value (mg/dL)	Relative Deviation (%)
1	85,0	85,77	-0,89
2	85,0	85,77	-1,04
3	85,0	85,77	-1,04

Precisión:

Within-run precision was determined with 3 lots and 10 replicates in an analytical run, between-run precision was determined with 3 lots and 3 replicates in an analytical run.

Samples (mg/dL)	Within-Run Precision		Between-Run Precision	
	SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
86.54	1.73	1.98	1.65	0.58

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES (THERAPEUTIC TARGET VALUES)

Suggested therapeutic target values for cardiovascular risk category estimated by the Doctor for adults over 20 years:

Risk	Therapeutical Target
Low	< 130 mg/dL
Intermediary	< 100 mg/dL
High	< 70 mg/dL
Very High	< 50 mg/dL

Conversion to the International System of Units (g/L):

LDL (mg/dL) x 0,026 = LDL (mmol/L)

NOTE: According to the Brazilian Consensus for the Normalization of Laboratory Determination of the Lipid Profile, reference values for LDL cholesterol were replaced by the values of therapeutic target, stratified according to the calculated cardiovascular risk.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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