Lactato

Lactate | Lactato
Ref. 10.018.00 / FD 38.073.00 / SD 39.070.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine lactate in plasma and cerebrospinal fluid samples. Diagnostic use only.



STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- The Dedicated and Semi Dedicated Reagents presentations have *onboard* stability of 30 days.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic - Trinder

The lactate in the sample, in the presence of oxygen, is oxidized by the enzyme lactate oxidase, generating pyruvate and hydrogen peroxide. The latter, in the presence of a phenolic reagent (TOOS), 4-aminoantipyrine and peroxidase produces a chromogen that is spectrophotometrically measured at 546 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: Plasma fluoride and cerebrospinal fluid.

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Plasma*	15 a 25 °C	8 hours
	4 to 8 °C	14 days
	-20 °C	1 month
Cerebrospinal Fluid**	4 to 8 °C	1 hour
	-20 °C	1 month

* Plasma samples should be separated within 15 minutes.

** Cerebrospinal fluid samples must be centrifuged.

PRODUCT DESCRIPTION

R 1

Pipes buffer ≥ 10 mmol/L; lactate oxidase ≥ 100 U/L; peroxidase ≥ 500 U/L; 4-aminoantipyrine ≥ 0,1 mmol/L; TOOS ≥ 0,5 mmol/L; cofactor; stabilizers; activators; detergent; preservative.

CAL

Phosphate buffer ≥ 20 mmol/L; bovine albumin ≥ 5% w/v; stabilizer; preservative and lactate at a concentration equivalent to 40 mg/dL. Traceable to NIST SRM 84L reference material and IC Assay NIST SRM.

CONTROL

Phosphate buffer ≥ 20 mmol/L; bovine albumin ≥ 5% w/v; stabilizer; preservative and lactate at the specified concentration on the vial label. Traceable to NIST SRM 84L reference material and IC Assay NIST SRM.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00
Normal Control - Quantinorm	13.003.00
Pathological Control - Quantialt	13.004.00

REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 546 nm (530 – 570 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Sample	Calibrator
CAL	---	---	10 µL
Sample	---	10 µL	---
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 5 minutes.

3. Measure the Sample's and Calibrator's absorbance against the Blank at 546 nm (530 nm – 570 nm).

B) CALCULATIONS

Lactate (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{CAL Absorbance}} \times \text{CAL Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{CAL Concentration (mg/dL)}}{\text{CAL Absorbance}}$

Lactate (mg/dL) = Sample Absorbance x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Plasma:

Lactate is an intermediate product from the metabolism of carbohydrates that predominantly accumulates in the blood when tissue oxygenation decreases. It originates mainly from the skeletal muscle, brain, skin, kidney marrow and erythrocytes. The concentration of lactate in the blood depends on the rate of production in these tissues and on the metabolic rate in the liver and kidneys. The liver uses approximately 65% (75 g/day) of total basal lactate, produced predominantly in gluconeogenesis. Anaerobic glycolysis markedly increases lactate in the blood, especially in prolonged exercise. The anoxia resulting from shock, pneumonia and congestive heart failure, Thiamine deficiency and diabetic ketoacidosis, among other causes, increase lactate and blood pyruvate. High concentrations of lactate in the blood result in lactic acidosis, which can occur in renal failure and leukemia.

CSF:

In cerebrospinal fluid (CSF), high lactate levels are used mainly in the differential diagnosis between bacterial and viral meningitis and are not limited to this purpose only. The increase in lactate can also occur due to any clinical condition that results in reduced oxygen flow to the tissues, such as destruction of the tissue within the Central Nervous System. In addition, lactate levels in CSF are often used to monitor serious head injuries.

PERFORMANCE CHARACTERISTICS

Operating range

Plasma	2.32 to 120.00 mg/dL
Cerebrospinal Fluid	2.07 a 80.00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity		
Sample Type	Detection Limit	Quantification Limit
Plasma	0.84 mg/dL	2.32 mg/dL
Cerebrospinal Fluid	-	1.87 mg/dL

Analytical Specificity			
Sample Type	Hemoglobin	Bilirubin	Triglycerides
Plasma	500 mg/dL	40 mg/dL	1000 mg/dL
Cerebrospinal Fluid	500 mg/dL	-	-

Concentrations of interfering substances up to the values presented above do not cause significant changes in the results. Any concentration of membrane stabilizer (Transfix*) interferes with the determination. For drugs, consult the recommended bibliographic reference (Young, 2000).

Sample Type	Accuracy	
	Number of samples	40 in duplicate
Plasma	Regression equation	y = 0.990x - 0.16
	Correlation coefficient (R)	0,9983
	Number of samples	20 in duplicate
Cerebrospinal Fluid	Regression equation	Y=0,9512x + 0,9141
	Correlation coefficient (R)	0,9976

Plasma: by applying the regression equation, the total systematic error estimated for decision levels of 50 mg/dL and 100 mg/dL was -1,32% and -1,16%, respectively.

Cerebrospinal Fluid: by applying the regression equation, the total systematic error estimated for decision levels of 8,50 mg/dL, 35,86 mg/dL and 75,00 mg/dL was 5,88%, -2,33% and -3,66%, respectively.

Sample Type	(mg/dL)	Replicates	Within-Run Precision		Total Precision	
			SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
Plasma	25.28	80	0.244	1.00	0.568	2.20
	35.25	80	0.247	0.70	0.524	1.50
	105.41	80	0.230	0.20	0.631	0.60
Cerebrospinal Fluid	16.20	20/40	0.398	2.46	0.328	2.02
	35.66	20/40	0.524	1.55	0.780	2.18
	84.29	20/40	1.308	1.55	1.452	1.73

%CV: Coefficient of Variation expressed as a percentage; SD: Standard Deviation

- RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**
- Use Protective Equipment in accordance with the Good Laboratory Practices.
 - Follow the Good Laboratory Practices' instructions to establish the quality of water.
 - Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
 - Avoid leaving reagents outside the specified storage conditions.
 - The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES		
	Plasma	CSF (Adults)
Venous	4.5 – 19.8 mg/dL	9.01 – 19.82 mg/dL
	0.5 – 2.2 mmol/L	
Arterial	4.5 – 14.4 mg/dL	1.00 – 2.20 mmol/L
	0.5 – 1.6 mmol/L	

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):
Lactate (mg / dL) x 0.111 = Lactate (mmol/L).

- WARNINGS AND PRECAUTIONS**
- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
 - The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotechnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotechnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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