Lipase VET

Lipase VET/ Lipasa VET
Ref. 90.050.00

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ENGLISH

INTENDED USE

Kit intended to determine the activity of Lipase in serum and plasma samples. Veterinary diagnostic use only.



STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: DGGR Colorimetric

Lipase catalyzes the cleavage of the 1-2-O-dilauril-rac-glycerol-3-glutaric-(6'-methylresorufin)-ester substrate in alkaline medium forming 1-2-O-dilauril-rac-glycerol and glutaric-6'-methylresorufin ester which spontaneously decomposes into glutaric acid and methylresorufin. The velocity of formation of methylresorufin is proportional to the catalytic activity of the lipase in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma heparin

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	3 weeks
	20 to 25 °C	7 days
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1

Bicine buffer ≥ 50 mmol/L; Colipase ≥ 1 mg / L; sodium taurodeoxycholate ≥ 5 mmol/L, sodium deoxycholate ≥ 1 mmol/L, activator, detergent, preservative.

X

R 2

Tartarato buffer ≥ 2 mmol/L; sodium taurodeoxycholate ≥ 5 mmol/L, Colipase $\geq 0,1$ mg / L; 1-2-O-dilauril-rac-glicerol-3-glutaric-(6'-methylresorufin)-éster ≥ 2 mmol/L; Stabilizers; activators, detergents, preservative.

X

Traceable to the specific absorbance of the chromogen Methylresorufin using calibrated pipettes and handheld spectrophotometer that provide absolute values.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the control below:

Calibrador – Autocal VET	REF	90.039.00
--------------------------	------------	-----------

Normal Control – Quantinorm VET	90.040.00
Pathological Control – Quantialt VET	90.041.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 580 (550 - 600) nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Adjust the temperature of the photometer to 37 °C, select the wavelength at 580nm and clear with purified water.
2. Pipette in the test tubes:

	Blank	Sample	Calibrator
Purified water	10 µL	-	-
Sample	-	10 µL	-
Calibrator	-	-	10 µL
R1	800 µL	800 µL	800 µL
R2	200 µL	200 µL	200 µL

3. Homogenize and insert in the thermostated cuvette. Set the stopwatch.
4. After 90 seconds, measure the initial absorbance (A1) and read it again in 180 seconds (A2).

B) CALCULATIONS

$$\text{Lipase (U/L)} = \frac{\Delta A \text{ Sample} - \Delta A \text{ Blank}}{\Delta A \text{ Calibrator} - \Delta A \text{ Blank}} \times \text{Concentration Calibrator (U / L)}$$

Where: $\Delta A = A2 - A1$

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Lipase is a low-molecular-weight protein, approximately 42 kDa, of high interest in the diagnosis of pancreatic diseases. Synthesized primarily in the pancreas, this enzyme plays a crucial role in the hydrolysis of triglycerides. The increase or decrease in lipase activity in animal serum may indicate specific clinical conditions, such as pancreatitis, pancreatic duct obstruction, pancreatic cancer, or cholelithiasis. On the other hand, decreases may be associated with exocrine pancreatic insufficiency. However, the specificity of assays may be affected by conditions such as renal insufficiency, enteritis, hepatopathies, bile duct carcinoma, and gastrointestinal tract lymphosarcoma. Therefore, the interpretation of lipase levels in the serum should consider the clinical context and be performed by healthcare professionals for an accurate diagnosis.

PERFORMANCE CHARACTERISTICS

Operating range
3,78 to 300 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
3,78 U/L

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	500 mg/dL	40 mg/dL	1000 mg/dL
Felines	300 mg/dL	10 mg/dL	1000 mg/dL
Equines	200 mg/dL	5 mg/dL	1000 mg/dL
Bovines	100 mg/dL	30 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy – Canines	
Number of Samples	20 in duplicate
Regression Equation	$y = 0,9918x - 0,65$
Correlation Coefficient (R)	0,9994
Accuracy – Felines	
Number of Samples	20 in duplicate
Regression Equation	$y = 0,9948x - 0,1687$
Correlation Coefficient (R)	0,9983
Accuracy – Equines	
Number of Samples	15 in duplicate
Regression Equation	$y = 0,9734x + 0,6583$
Correlation Coefficient (R)	0,9947
Accuracy – Bovines	
Number of Samples	10 in duplicate

Regression Equation	y = 1.0769x - 0.0603
Correlation Coefficient (R)	0.9885

Precision:

Determined with two runs.

Samples (U/L)	Within-Run Precision		Samples (U/L)	Total Precision	
	SD (U/L)	CV (%)		SD (U/L)	CV (%)
4.91	0.12	2.4	5.03	0.22	4.5
133.50	2.06	1.5	133.25	2.42	1.8

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Canines	13 – 200 U/L
Felines	0 – 83 U/L
Equines	7 – 16 U/L
Bovines	5 – 13 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

Lipase (U/L) x 0,0167 = Lipase (µkat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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