## Lisante HemaBio V

Lyse HemaBio V | Lisante HemaBio V  
Ref. 14.033.00

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### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

### INTENDED USE

Lysing solution for hematological analyzers. Use for *in vitro* diagnostics only.

### STORAGE AND HANDLING

- Ready-to-use reagents.
- Store at 2 to 45 ° C. Do not freeze.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Once opened, the product is stable for 4 months, as long as the recommended storage conditions (2 to 45°C) are followed.
- Do not use reagents whose date has expired.

### WORKING PRINCIPLE

The hematological lyser promotes rapid lysis of red blood cells (erythrocytes), releasing hemoglobin and changing the apparent volume of white blood cells. When cells pass through the opening of the electrode present in hematological analyzers, they generate electrical pulses due to the difference in conductivity between the diluent and the cells. These pulses, as established by Coulter's principle, are proportional to the quantity and size of blood cells, allowing for cell sorting and counting.

### PRODUCT DESCRIPTION

Quaternary solution of ammonium salt.

### TEST PROCEDURE

Open the outer packaging or plastic lid and sealed membrane and insert the hematology analyzer conduit into the reagent. For more operation details, please refer to the hematology analyzer user manual. Leave the reagent at room temperature for at least 24 hours before use.

### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Avoid contact with skin and eyes as quaternary ammonium salts may cause irreversible effects to these tissues.
- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

### ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone +55 35 3214 4646.
- Dispose reactions waste in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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