	
<h2>Magnésio Mono VET</h2>	Responsável Técnico: Dr. Gilson Sérgio Pizzo CRF MG – 5310
Magnésio Mono VET / Magnésio Mono VET Ref. 90.026.00	

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ENGLISH

INTENDED USE
 Kit intended to determine magnesium in serum and plasma samples. Veterinary Diagnostic use only.

STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Sulfonated Magon
 Magnesium forms with sulfonated Magon (xylydid blue I), in alkaline medium, a colored complex that can be spectrophotometrically measured at 505 nm. Color intensity is proportional to the magnesium concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum and heparin plasma.
Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

	Temperature	Sample Stability
Serum and plasma	4 to 8 °C -20 °C	7 days 6 months

PRODUCT DESCRIPTION

R 1 Tris buffer ≥ 100 mmol/L; sulphonated magon ≥ 0.1 mmol/L; activators; preservative; stabilizers. ✘

STD Magnesium ions in concentration equivalent to 2,0 mg/dL; preservative. Traceable to the reference material NIST 3131a. ✘

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Autocal VET	REF	90.039.00
Normal Control – Quantinorm VET		90.040.00
Pathological Control – Quantialt VET		90.041.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 505 nm (500 – 510 nm).
- Test tubes, glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 15 to 30 °C for 2 minutes.

3. Measure the Sample's and Standard's absorbance against the blank at 505 nm (500 – 510 nm). The final reaction is stable for 15 minutes.

B) CALCULATIONS

Calculation of serum and plasma
 Magnesium (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$
 Magnesium (mg/dL) = Sample Absorbance x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaltda.ind.br

C) INTERPRETATION

Magnesium is an essential mineral for the proper functioning of the animal's body. It acts in more than 300 enzymatic processes as a cofactor, allosteric activator and/or in the formation of enzyme substrates, such as MgATP. Ingested magnesium is absorbed mainly in the small intestine. Its regulation occurs in the kidneys, in the ascending branch of Henle. Hypomagnesemia can be observed in several clinical conditions, such as: anorexia, gastrointestinal disorders, acute pancreatitis, cholestasis, glomerulonephritis, prolonged intravenous fluid therapy, diabetic ketoacidosis, hyperthyroidism, sepsis, total parenteral nutrition. Hypermagnesemia can occur in some situations, for example: excessive use of antacids, enemas and parenteral fluids containing magnesium, rhabdomyolysis, acute and chronic renal failure.

PERFORMANCE CHARACTERISTICS

Operating range
0.40 to 4.50 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.07 mg/dL	0.40 mg/dL

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	300 mg/dL	30 mg/dL	1000 mg/dL
Felines	100 mg/dL	40 mg/dL	750 mg/dL
Equines	300 mg/dL	10 mg/dL	1000 mg/dL
Bovines	300 mg/dL	10 mg/dL	750 mg/dL

Concentrations of interfering substances up to the values presented above do not cause significant changes in the results. Any concentration of membrane stabilizer (Transfix®) interferes with the determination. For drugs, consult the recommended bibliographic reference (Young, 2000).

Accuracy		
Control	Control 1	Control 2
Interval (mg/dL)	1.98 – 2.77 (target: 2.375)	3.9– 4.33 (target 4.12)
Average (mg/dL)	2.35	4.14
Replicates	10	10
%DIF	-0.8	-4.4

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
1.32	80	0.008	0.60	0.009	0.70
3.45	80	0.034	1.00	0.039	1.10
4.30	80	0.049	1.10	0.056	1.30

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results

- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs of disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

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These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Magnesium (mg/dL) x 0,41 = Magnesium (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br

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