Magnésio Mono

Magnesium Mono / Magnésio Mono
Ref. 12.007.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine magnesium in serum, plasma and urine samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Sulfonated Magon

Magnesium forms with sulfonated Magon (xylylid blue I), in alkaline medium, a colored complex that can be spectrophotometrically measured at 505 nm. Color intensity is proportional to the magnesium concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum, heparin plasma and urine.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preparation:

	Temperature	Sample Stability
Serum and Plasma	4 to 8 °C	7 days
	-20 °C	6 months
24 Hours Urine	4 to 8 °C pH < 2,0	3 days

Preparation:

Urine: work with samples collected in the period of 24 hours in a flask with 20 to 30 mL of a 50% HCl solution. Homogenize the urine, measure its volume and adjust the pH to approximately 1,0. Perform a 1:5 dilution of an urine aliquot with purified water and perform the assay. Multiply the result obtained by 5. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

PRODUCT DESCRIPTION

R 1	Tris buffer ≥ 100 mmol/L; sulphonated magon ≥ 0.1 mmol/L; activators; preservative; stabilizers. 
STD	Magnesium ions in concentration equivalent to 2,0 mg/dL; preservative. Traceable to the reference material NIST 3131a. 

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Autocal H	13.002.00	90.00
Normal Control – Quantinorm	13.003.00	90.00
Pathological Control – Quantialt	13.004.00	90.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 505 nm (500 – 510 nm).
- Test tubes, glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 15 to 30 °C for 2 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 505 nm (500 – 510 nm). The final reaction is stable for 15 minutes.

B) CALCULATIONS

Calculation of serum and Plasma

Magnesium (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Magnesium (mg/dL) = Sample Absorbance x Calibration Factor

Calculations for urine samples:

Magnesium urine (mg/dL) = $\frac{\text{Sample absorbance}}{\text{Standard Absorbance}} \times \text{STD concentration (mg/dL)} \times 5$

*The urine sample must be diluted in a 1:5 proportion. To determine the magnesium concentration in the urine sample (mg/dL), multiply the obtained result by 5.

Magnesium (mg/24 horas) = $\frac{\text{Magnesium (mg/dL)} \times \text{Urinary Volume in 24 hours (mL)}}{100}$

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Approximately 55% of the human Mg²⁺ cation is found in the skeleton, another part inside tissue cells and a small proportion in the extracellular fluid (about 1%), being the fourth most abundant cation in the body. Magnesium ion acts in more than 300 enzymatic processes as cofactor, allosteric activator and / or in the formation of enzymatic substrates, such as MgATP. Ingested magnesium is absorbed primarily in the small intestine. Its regulation occurs in the kidneys, in the ascending branch of Henle. Hypomagnesemia can be observed in several clinical conditions, such as: alcoholism, acute pancreatitis, malabsorption syndromes, metabolic acidosis, renal diseases and protein-calorie malnutrition. Hypermagnesemia may occur in some situations, for example: excessive use of antacids, enemas and parenteral fluids containing magnesium, rhabdomyolysis, acute and chronic renal failure. The dosage of urinary magnesium is mainly used to evaluate hydroelectrolytic.

PERFORMANCE CHARACTERISTICS

Operating range
0,40 to 4,50 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.07 mg/dL	0.40 mg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
500 mg/dL	40 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1.019x - 0.034$
Correlation Coefficient (R)	0.9940

By applying the regression equation, the total systematic error estimated for decision levels of 2.0 mg/dL and 3.0 mg/dL was 0.20% and 0.77%, respectively.

Precisión:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
1.32	80	0.008	0.60	0.009	0.70
3.45	80	0.034	1.00	0.039	1.10
4.30	80	0.049	1.10	0.056	1.30

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Serum and Plasma	1,6 - 2,6 mg/dL
Urine 24h	72,9 – 121,5 mg/24 hours

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Magnesium (mg/dL) x 0,41 = Magnesium (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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