Microalbuminúria

Microalbumin / Microalbuminúria
Ref. 20.014.00 (Microalbuminúria)
Ref. 20.014.00 (Microalbuminúria - C)

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the quantitative determination of albumin in the urine. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagents, control and calibrator.

WORKING PRINCIPLE

Method: Immunoturbidimetric

Latex particles coated with anti-human albumin antibodies are agglutinated by the albumin present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of albumin.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: Urine

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Urine	2 to 8 °C	7 days*

*When sodium azide 1 g / L is added.

PRODUCT DESCRIPTION

R DIL	100 mmol / L glycine buffer pH 10.0, preservative.		
R LAX	Latex particles coated with goat anti-human albumin IgG antibodies pH 8.2, preservative.		
CAL	Human Albumin, preservative. The albumin concentration is traceable to the International Reference Material CRM 470 / RPPHS of BCR.		
CONTROL	Solution of human Albumin, fosfat buffer, stabilize, preservative.		

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Microalbumin CAL (Kit Included)	REF	20.014.00
Microalbumin CONTROL (Kit in Microalbum - C)		20.014.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 540 nm (520-560 nm).

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

REACTION PARAMETERS	
Type of reaction	Fixed time
Wave length	540 nm (520-560 nm)
Temperature	37 °C
Sample volume/CAL*	7,0 µL
RDIL VOLUME*	800 µL
RLAX VOLUME *	200 µL
Reading 1 (A1)	Immediately
Reading 2 (A2)	2 min.

* Volumes can be modified as long as the established ratio is maintained.

B) CALCULATIONS

Albumin (mg/L) = $\frac{(A2 - A1) \text{ of Sample} \times \text{CAL Concentration (mg/L)}}{(A2 - A1) \text{ of CAL}}$

Albumin (mg/24 hours) = Albumin (mg/L) x 24-hour urine volume (L)

Albumin Rate/Creatinine* (mg/g) = $\frac{\text{Albumin (mg/L)}}{\text{Urinary creatinine (g/L)}}$

* Albumin and creatinine should be determined in the same sample.

C) INTERPRETATION

Microalbumin is an increase in urinary albumin excretion as a consequence of changes in glomerular permeability. Increased urinary albumin excretion precedes and is highly predictive of diabetic nephropathy, end-stage renal disease and proliferative retinopathy in type 1 diabetes; it is a predictor of progressive renal disease in patients with type 2 diabetes; and identifies a group of non-diabetic individuals at high risk for coronary artery disease.

PERFORMANCE CHARACTERISTICS

Operating range
5,1 a 163 mg/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Limit of detection	Limit of quantification
1.07 mg/L	5.10 mg/L

Analytical Specificity			
Hemoglobin	Bilirubin	Urea	Creatinine
100 mg/dL	10 mg/dL	150 mg/dL	300 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	39 in duplicate
Regression Equation	$y = 1.07x + 0.587$
Correlation Coefficient (R)	0.9995

By applying the regression equation, the total systematic error estimated for decision levels of 60 mg/L and 130 mg/L was 1.08% and 0.55%, respectively.

Prozone High Dose Effect:

There is no prozone effect up to 402 mg/L of Albumina

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/L)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/L)	CV (%)	SD (mg/L)	CV (%)
22,3	80	0,497	2,2	0,638	2,9
51,9	80	1,304	2,5	1,335	2,6
113,4	80	2,279	2,0	2,735	2,4

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

	Albumin (mg/24hs)	Albumin/Creatinine (mg/g)
Normal to slightly increased	< 30	< 30
Moderately increased	30 a 300	30 a 300
Severely increased	> 300	> 300

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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- BAR, J. et al., Microalbuminuria: Prognostic and Therapeutic Implications in Diabetic and Hypertensive Pregnancy. Diabetic Medicine. v.12, p. 649-656, 1995.
- CAMBIASSO, C.L. et al., Immunoassay of Low Concentrations of Albumin in Urine by Latex Particle Counting. Clin. Chem. v. 34, p. 416-418, 1988.
- Kidney Disease: Improving Global Outcomes (KDIGO) CKD Work Group. KDIGO clinical practice guideline for the evaluation and management of chronic kidney disease. Kidney Int. Suppl. 2013, 31-150.
- RASMUSSEN, B.F; et al., Microalbuminuria: An Important Diagnostic Tool. J Diab Comp. v. 8, p. 137-145, 1994.
- MEDCALF, E.A; et al., Rapid, Robust Method for Measuring Low Concentrations of Albumin in Urine. Clin. Chem. v. 36, p. 446-449, 1990.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27, p.493-501, 1981.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

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