## PCR Turbilátex

CRP Turbilátex / PCR Turbilátex  
Ref. 20.015.00 / FD 38.054.00 / SD 39.063.00

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#### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

#### INTENDED USE

Kit for the determination of C-reactive protein (CRP) in serum. *In vitro* diagnostic use.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 28 days.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagents and calibrator

#### WORKING PRINCIPLE

**Method:** Immunoturbidimetric  
Latex particles coated with anti-human CRP antibodies are agglutinated by the CRP present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of CRP.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum	2 to 8 °C	7 days
	-20 °C	3 months

#### PRODUCT DESCRIPTION

<b>R DIL</b>	TRIS buffer, 20 mmol / L, pH 8.2, preservative.		
<b>R LAX</b>	Latex particles coated with goat anti-PCR IgG pH 7.3; preservative.		
<b>CAL</b>	Human serum, preservative. The CRP value is traceable to Reference Material CRM 470/RPPHA (BCR).		

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

CRP Turbilátex CAL (Kit Included)	20.015.00
Rheumatic Control Level 1	21.007.00
Rheumatic Control Level 2	21.008.00

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 540 nm (520 - 560 nm).

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

**Automação:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

REACTION PARAMETERS	
<b>Type of reaction</b>	Fixed time
<b>Wave length</b>	540 nm
<b>Temperature</b>	37 °C
<b>Sample volume/CAL*</b>	2 µL
<b>RDIL*</b>	320 µL
<b>RLAX*</b>	80 µL
<b>Reading 1 (A1)</b>	Immediately
<b>Reading 2 (A2)</b>	2 min.

\* Sample volumes/CAL, RDIL and RLAX be modified as long as the established ratio is maintained.

##### B) CALCULATIONS

CRP(mg/L) =  $(A_2 - A_1)$  of Sample x Concentration of CAL (mg/L) /  $(A_2 - A_1)$  of CAL

#### C) INTERPRETATION

C-Reactive Protein is an acute phase protein that can increase significantly in most bacterial and viral infectious processes, inflammation and malignant neoplasms. Its elevation occurs rapidly and can reach values of 300 mg / L in 12 to 24 hours.

#### PERFORMANCE CHARACTERISTICS

Operating range
3,12 to 123 mg/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0,38 mg/L	3,12 mg/L

Analytical Specificity			
Hemoglobin	Bilirubin	Triglycerides	Rheumatoid Factor
500 mg/dL	30 mg/dL	1000 mg/dL	600 UI/mL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1.006x - 0.048$
Correlation Coefficient (R)	0.9996

By applying the regression equation, the total systematic error estimated for decision levels of 25 mg/L and 150 mg/L was 0.41% and 0.57%, respectively.

#### Prozone High Dose Effect:

There is no prozone effect up to 297 mg/L of CRP.

#### Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/L)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/L)	CV (%)	SD (mg/L)	CV (%)
11,1	80	0,375	3,4	0,375	3,4
54,5	80	1,787	3,3	1,991	3,7
73,2	80	1,470	2,0	1,836	2,5

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

#### REFERENCE RANGES

Adults	Up to 5 mg/L
--------	--------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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- ADAM, A. *et al.*, Analytical Evaluation of a New Latex Agglutination Test for Quantitative Determination of C-Reactive Protein by Laser Nephelometry. *J. Clin. Chem. Clin. Biochem.* v. 23, p. 787-789, 1985.
- HOKAMA, Y.; NAKAMURA, R.M. C- Reactive Protein: Current Status and Future Perspectives. *Journal of Clinical Laboratory Analysis.* v. 1, p. 15-27, 1987.
- WESTGARD, J. O. *et al.* A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27, p.493-501, 1981
- WINKLES, J. *et al.*, Enhanced-Latex-Agglutination Assay for C-Reactive Protein in Serum, with Use a Centrifugal Analyzer. v. 33, p. 685-689, 1987.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

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