Sensitivity	
Quantification Limit	
0,40 mmol/L	

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	500 mg/dL	10 mg/dL	250 mg/dL
Felines	200 mg/dL	18 mg/dL	500 mg/dL
Equines	300 mg/dL	10 mg/dL	250 mg/dL
Bovines	100 mg/dL	10 mg/dL	250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy				
	Caninos	Felinos	Equinos	Bovinos
Number of Samples	11 in duplicate	10 in duplicate	8 in duplicate	8 in duplicate
Regression Equation	y = 1.0232x – 0.1045	y = 1.0372x – 0.1896	y = 1.0102x – 0.1019	y = 0,9619x+ 0,2147
Correlation Coefficient (R)	0,9755	0,9814	0,9964	0,9771

Precision:

Samples (mmol/L)	Repetition s	Within-Run Precision		Samples (mmol/L)	Repetition s	Total Precision	
		SD (mg/dL)	CV (%)			SD (mg/dL)	CV (%)
4,40	20	0,21	4,8	4,34	40	0,19	4,3
6,99	20	0,18	2,5	6,95	40	0,27	3,9

Determined with two runs in duplicate.

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the requirements by means of periodic checks.

REFERENCE RANGES

Canine	4,37 – 5,35 mmol/L
Feline	4,0 – 4,5 mmol/L
Equine	2,4 – 4,4 mmol/L
Bovine	3,9 – 5,8 mmol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Potassium (mmol/L) x 1.0 = Potassium (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechind.com.br.

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