---

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Sensitive Protein | Proteína Sensível  
Ref. 10.016.00 (KIT) / 10.016.00 (KIT COM CONTROLE)

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**ENGLISH**

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

**INTENDED USE**

Kit intended for the quantitative determination of protein in urine and cerebrospinal fluid (CSF) samples. *In vitro* diagnostic use.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**

**Method:** Pyrogallol Red

The protein reacts with pyrogallol red and sodium molybdate, in acidic medium, forming a colored complex that can be spectrophotometrically measured at 600 nm. The intensity of the color is proportional to the concentration of protein in the sample.

**SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY**

**Sample Type:** urine and cerebrospinal fluid (CSF).

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Urine	4 to 8 °C	7 days
	-20 °C	1 month
Cerebrospinal Fluid	4 a 8 °C	6 days
	-20 °C	1 year

**Preparation:**

	Temperature	Stability Period
Urine	4 to 8 °C	7 days
	-20 °C	1 month
Cerebrospinal Fluid	4 a 8 °C	6 days
	-20 °C	1 year

---

**URINE:** homogenize the sample and separate 20 mL aliquot. Centrifuge for 10 minutes at 3000 rpm and utilize the supernatant to proceed with the assay.

**Cerebrospinal fluid:** centrifuge for 10 minutes at 3000 rpm and use the supernatant in the assay.

**PRODUCT DESCRIPTION**

<b>R 1</b>	Succinate buffer > 10 mmol / L; sodium molybdate > 10 µmol / L; pyrogallol red > 10 µmol / L; preservative; detergent; stabilizers. Pipes buffer; preservatives; stabilizer; Bovine albumin at a concentration corresponding to 50,0 mg / dL protein. Traceable to NIST reference material 927d.
<b>STD</b>	
<b>CONTROL</b>	Phosphate buffer; stabilizers; detergent; preservatives and bovine albumin. (INCLUDED IN THE KIT WITH CONTROL)

---

**QUALITY CONTROL**

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the control below:

SENSIBLE PROTEIN CONTROL (INCLUDED IN THE KIT WITH CONTROL)	REF	10.016.00
---	-----	-----------

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 600 nm (580 – 620 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Pipette in the test tubes:

	Blank	Sample	STD
STD	---	---	50 µL
Sample	---	50 µL	---
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate for 10 minutes at 37 °C. Remove it from the water bath and let it rest in temperature of 15 to 30 °C before reading.

3. Measure the Standard's and Sample's absorbance against the blank at 600 nm (580 – 620 nm). The final reaction is stable for 15 minutes.

---

**B) CALCULATIONS**

Protein (mg/dL) =  $\frac{\text{Sample Absorbance} - \text{Blank}}{\text{STD Absorbance} - \text{Blank}} \times \text{STD Concentration (mg/dL)}$

**Calculations with the Calibration Factor:**

Calibration Factor =  $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance} - \text{Blank}}$

Protein (mg/dL) = (Sample Absorbance - Blank) x Calibration Factor

**Calculations for urine samples:**

Protein (mg/24 hours) =  $\frac{\text{Protein (mg/dL)} \times \text{Urinary Volume (mL)}}{100}$

**Automação:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

---

**C) INTERPRETATION**

**Urine:**

Proteinuria is defined as the presence of excess protein in the urine, which occurs in many types of kidney disease, after intense exercise or dehydration. The causes may be due to increased glomerular permeability, altered tubular reabsorption, elevated low-molecular-weight plasma proteins, abnormal secretion of proteins in the urinary tract and post-renal bleeding.

**Cerebrospinal fluid:**

Examination of CSF total proteins is mainly used to detect increased permeability of the blood-brain barrier or increased intrathecal production of immunoglobulins. The permeability of the blood-brain barrier to plasma proteins is increased by intracranial pressure due to brain tumor, intracerebral hemorrhage, traumatic injury, inflammation due to bacterial or viral meningitis, encephalitis or poliomyelitis.

**PERFORMANCE CHARACTERISTICS**

Sensitivity		
Sample type	Limit of Detection	Limit of Quantification
Urine	0.50 mg/dL	4.80 mg/dL
Cerebrospinal fluid	-	14.39 mg/dL

Operating range	
Urina	4.80 a 97.00 mg/dL
Cerebrospinal fluid	14.28 a 68.55 mg/dL

---

**BIOTÉCNICA IND.COM. LTDA - CNPJ: 025340690001/20.** Avenida Washington Ribeiro, 200 - Industrial Miguel de Luca CEP 37072-030 Varginha-MG BRASIL Tel/fax: +55 35 3214 4646 [www.biotechnica.ind.br](http://www.biotechnica.ind.br)

Revisão/Revision/Revisión: 04 – 04/04/2024