CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices’ instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES	
Prothrombin time	11 to 15 seconds
RNI - healthy people	0.9 to 1.1
RNI - patients receiving oral anticoagulant therapy	2.5 to 3.5

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- GUYTON, A. C.; HALL, J. E. *Fisiologia humana e mecanismos das doenças*. 6ª ed. Rio de Janeiro: Guanabara Koogan. 1998. 580p.
- BURTIS, Carl A.; ASHWOOD, Edward R. **Tietz: Fundamentos de Química Clínica**. 4. ed. Philadelphia: W. B. Saunders, 1996. 836 p.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. **Tietz: Fundamentos de Química Clínica**. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.
- CLSI. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline-Second Edition; *Approved Guideline. CLSI document H47-A2*. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- HORSTI J.; UPPA H.; VILPO J.A. Poor Agreement among Prothrombin Time International Normalized Ratio Methods: Comparison of Seven Commercial Reagents. **Clin. Chem.** V. 51 N° 3, p.553-560, 2005.
- POLLER, L. International Normalized Ratios (INR): the first 20 years. **Journal Of Thrombosis And Haemostasis**. Manchester, p. 849-860, 2004.
- HARRIS, Neil. The International Normalized Ratio: How well do we understand this measurement? **American Association for Clinical Chemistry: AACC**, 04 set. 2012.

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