**QUANTIALT**

QUANTIALT / QUANTIALT  
Ref. 13.004.00

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**ENGLISH**

**INTENDED USE**

Kit intended for the internal quality control of quantitative biochemical assays. Use in *in vitro* diagnostic.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- After reconstituted, the analytes are stable for: Bilirubins and Alkaline Phosphatase: 48 hours when stored at 2 to 8°C, in the original flask. Other analytes: 5 days when stored at 2 to 8°C, in the original flask, and 15 days when stores at -20 °C, in microtubes.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**

Human lyophilized serum used in the internal quality control of biochemical assays.

**PREPARATION AND HANDLING OF THE PRODUCT**

Reconstitute with 5 mL of purified water and let stand for 30 minutes. Gently homogenize by inversion, avoiding foaming. Until the complete dissolution of control.

**PRODUCT DESCRIPTION**

**CONTROL** Human serum plus tissue extracts of human/or animal origin. 

The control shows negative results for the HBs antigen and for the HIV and HCV antibodies. However, they should be treated as potentially infectious. The control values are traceable to reference methods and reference materials from the National Institute of Standard and Technology (NIST) and European Reference Materials (ERM).

**NECESSARY EQUIPMENT FOR TESTING**

- Glass and/or automatic pipettes.
- Purified water.

**TEST PROCEDURE**

The control should be used in the same way as patient samples, according to the procedure established for the assays and in accordance with the equipment protocol.

**QUALITY CONTROL**

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision

level and another control with value in another range of clinical significance. For internal control of laboratory quality it is recommended to use the controls below:

QUANTINORM Control	<b>REF</b>	13.003.00
QUANTIALT Control		13.004.00

**RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix controls from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the product if it displays any signs in disagreement with the ones specified in its MSDS.
- Avoid leaving controls outside the specified storage conditions.

**WARNINGS AND PRECAUTIONS**

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646.
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper lace for potentially infectious material.

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnicatda.com.br](mailto:sac@biotecnicatda.com.br)

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- CEMBROWSKI, G. Thoughts on quality-control systems: A laboratorians perspective. Clin. Chem. v.43, p.886-892, 1997.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981. Approved Guideline. **CLSI document H47-A2**. Wayne, PA: Clinical and Laboratory Standards Institute; 2008

Lote / Lot / Lote: 10M080		Vencimento / Due Date / Validez: 03/2026		
Analito / Parameter / Analito	Metodologia / Method / Método	Média / Average / Media	Intervalo / Range / Intervalo	Unidade / Unity / Unid
Ácido Úrico / Uric acid / Ácido Úrico	Enzimático Colorimétrico / Enzimatic Colorimetric / Enzimático Colorimétrico	8,57	7,11 – 10,0	mg/dL
Albumina / Albumin / Albúmina	Verde de Bromocresol / Bromocresol Green / Verde de Bromocresol	4,06	3,65 – 4,47	g/dL
Bilirrubina Direta / Direct Bilirubin / Bilirrubina Directa	DCA - Dicloroanilina Diazotada / Diazotized Dichloroaniline / Dicloroanilina Diazotada	3,74	2,99 – 4,49	mg/dL
	Dimetilsulfóxido / dimethylsulfoxide / dimetilsulfóxido	3,87	3,10 – 4,64	mg/dL
Bilirrubina Total / Total Bilirubin / Bilirrubina Total	DCA - Dicloroanilina Diazotada / Diazotized Dichloroaniline / Dicloroanilina Diazotada	5,40	4,32 – 6,48	mg/dL
	Dimetilsulfóxido / dimethylsulfoxide / dimetilsulfóxido	5,21	4,17 – 6,25	mg/dL
Creatinina / Creatinine / Creatinina	Picrato Alcalino / Alkaline Picrate / Picrato Alcalino	4,66	3,96 – 5,36	mg/dL
	Picrato Alcalino Monorreag. / Alkaline Picrate Monoreagent / Picrato Alcalino Reactivo único	4,71	4,00 – 5,42	mg/dL
	Enzimático / Enzimatic / Enzimático	4,99	4,24 – 5,74	mg/dL
Frutosamina / Fructosamine / Fructosamina	Azul de Nitrotetrazol / Nitrotetrazol Blue / Azul de nitro-tetrazolio <b>Sistemas Beckman AU</b>	393	314 – 472	µmol/L
	Azul de Nitrotetrazol / Nitrotetrazol Blue / Azul de nitro-tetrazolio <b>Demais Analisadores</b>	376	301 – 451	µmol/L
Glicose / Glucose / Glucosa	Enzimático Colorimétrico / Enzimatic Colorimetric / Enzimático Colorimétrico	288	259 – 317	mg/dL
Glicose Plus / Glucose Plus / Glucosa Plus	Enzimático GOD / Enzimatic GOD / Enzimático GOD	287	258 – 316	mg/dL
Lactato / Lactate / Lactato	Enzimático Colorimétrico / Enzimatic Colorimetric / Enzimático Colorimétrico	24,9	19,9 – 29,9	mg/dL
Proteínas Totais / Total Protein / Proteína Total	Colorimétrico / Colorimetric / Colorimétrico	6,73	6,06 – 7,40	g/dL
Ureia / Urea / Urea	Enzimático Colorimétrico / Enzimatic Colorimetric / Enzimático Colorimétrico	94,5	77,5 – 112	mg/dL
	Enzimático UV / Enzimatic UV / Enzimático UV	93,0	76,3 – 110	mg/dL

Colesterol / Cholesterol / Colesterol	Enzimático Colorimétrico / Enzimatic Colorimetric / Enzimático Colorimétrico	291	262 – 320	mg/dL
HDL Colesterol Direto / HDL-Cholesterol Direct / HDL Colesterol Directo	Acelerador Detergente Seletivo / Selective detergent accelerator / Acelerador Detergente seletivo	121	84,7 – 157	mg/dL
	Direto – Inibição seletiva / Direct - Selective inhibition / Directa - Inhibición selectiva	99,8	69,9 – 130	mg/dL
LDL Colesterol / LDL Cholesterol / LDL Colesterol	Surfactante Seletivo / Selective Surfactant / Surfactante seletivo	118	90,9 – 145	mg/dL
	Plus: Direto / Direct / Directo	157	121 – 193	mg/dL
Triglicérides / Triglycerides / Triglicéridos	Enzimático Colorimétrico / Enzimatic Colorimetric / Enzimático Colorimétrico	230	173 – 288	mg/dL
Alfa Amilase / Alpha Amylase / Alfa Amilasa	Cinético (CNP) / Kinetic (CNP) / Cinético (CNP)	499	349 – 649	U/L
CK-NAC / CK-NAC / CK-NAC	Cinético UV / Kinetic UV / Cinético UV	305	214 – 397	U/L
Colinesterase / Cholinesterase / Colinesterasa	Butirilcolina Cinético / Butyrylthiocholine Kinetic / Butirilcolina Cinética	4665	3732 – 5598	U/L
Desidrogenase Láctica / Lactate Dehydrogenase / Lactato Desidrogenasa	Cinético UV (P → L) / Kinetic UV (P → L) / Cinético UV (P → L)	783	626 – 940	U/L
Fosfatase Alcalina / Alkaline phosphatase / Fosfatasa Alcalina	p - nitrofenilfosfato cinético / Kinetic p-nitrophenylphosphate / p - nitrofenilfosfato cinético (DGKC)	336	235 – 437	U/L
	Cinético / Kinetic / Cinético (IFCC)	247	173 – 321	U/L
Gama GT / Gamma-glutamyl transferase / Gama GT	Cinético / Kinetic / Cinético	149	104 – 194	U/L
Lipase / Lipase / Lipasa	DGGR Colorimétrico / DGGR Colorimetric / Colorimétrico DGGR	118	94,4 – 142	U/L
TGO/AST / AST/GOT / GOT/AST	Cinético UV / Kinetic UV / Cinético UV	254	203 – 305	U/L
TGP/ALT / ALT/GPT / GPT/ALT	Cinético UV / Kinetic UV / Cinético UV	137	110 – 164	U/L
Cálcio / Calcium / Calcio	Arsenazo III / ASX / Arsenazo III	12,9	11,0 – 14,8	mg/dL
Cloreto / Chloride / Cloruros	Tiocinato / Thiocyanate / Toocianato	107	93,1 – 121	mEq/L
Cobre / Copper / Cobre	Colorimétrico / Colorimetric / Colorimétrico	90,1	72,1 – 108	µg/dL
IBC / IBC / UIBC	Ferrozine / Ferrozine / Ferrozina	215	172 – 258	µg/dL
Fósforo / Phosphorus / Fósforo	UV / UV / UV	7,47	6,72 – 8,22	mg/dL

Ferro / Iron / Hierro	Cromazurol / CRX / Cromazurol	105	84,0 – 126	µg/dL
	Ferrozine / Ferrozine / Ferrozine	235	188 – 282	µg/dL
Magnésio / Magnesium / Magnesio	Magon Sulfonado / Magon Sulfonated / Magon Sulfonado	3,64	2,73 – 4,55	mg/dL
Potássio / Potassium / Potasio	Enzimático UV / UV Enzimatic / Enzimático UV	7,96	7,16 – 8,76	mmol/L
Sódio / Sodium / Sodio	Enzimático / Enzimatic / Enzimático UV	170	153 – 187	mmol/L
Zinco / Zinc / Zinc	Colorimétrico / Colorimetric / Colorimétrico	33,8	27,0 – 40,6	µg/dL

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The values reported for the control serve as a reference. Each laboratory should establish the range of variation in its own analytical system.

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