QUANTIALT VET

QUANTIALT VET / QUANTIALT VET
Ref. 90.041.00

FINALIDADE

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ENGLISH

INTENDED USE

Kit intended for the internal quality control of quantitative biochemical assays. Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- After reconstituted, the analytes are stable for: Alkaline Phosphatase: 48 hours when stored at 2 to 8°C, in the original flask. Other analytes: 5 days when stored at 2 to 8°C, in the original flask, and 15 days when stores at -20 °C, in microtubes. Once defrosted, the analytes are stable for 24 hours.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Human lyophilized serum used in the internal quality control of biochemical assays.

PREPARATION AND HANDLING OF THE PRODUCT

Reconstitute with 5.0 mL of purified water. Homogenize gently by inversion, avoiding foaming. Wait 30 minutes until complete dissolution of the product.

PRODUCT DESCRIPTION

CONTROL Human serum plus tissue extracts of human/animal origin. 

The control shows negative results for the HBs antigen and for the HIV and HCV antibodies. However, they should be treated as potentially infectious.

NECESSARY EQUIPMENT FOR TESTING

- Glass and/or automatic pipettes.
- Purified water.

TEST PROCEDURE

The control should be used in the same way as patient samples, according to the procedure established for the assays and in accordance with the equipment protocol.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.

- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix controls from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the product if it displays any signs in disagreement with the ones specified in its MSDS.
- Avoid leaving controls outside the specified storage conditions.

WARNINGS AND PRECAUTIONS

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646.
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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- CEMBROWSKI, G. Thoughts on quality-control systems: A laboratorians perspective. Clin. Chem. v.43, p.886-892, 1997.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.
- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. **Clinical biochemistry of domestic animals**. 5th. ed. San Diego: Academic Press, c1997 932p.
- THRALL, Mary Anna. **Hematologia e bioquímica clínica veterinária**. São Paulo: Roca, 2007. 582 p.
- KERR, Morag G. Exames laboratoriais em medicina veterinária: **bioquímica clínica e hematologia**. 2. ed. São Paulo: Roca, 2003. 436 p.
- DUNCAN, J. Robert; PRASSE, Keith W. **Patologia clínica veterinária**. Rio de Janeiro: Guanabara Koogan, 1982. 217p.
- BUSH, B. M. **Manual del laboratorio veterinario de analisis clinicos**. Zaragoza: Acribia, 1982. 467p.
- Stockham, Steven L., Scott, Michael A. **Fundamentals of veterinary clinical pathology**. Blackwell Publishing, 2008.
- WILLARD, Michael D; TVEDTEN, Harold; TURNWALD, Grant H. **Small animal clinical diagnosis by laboratory methods**. 2nd ed. Philadelphia: W.B. Saunders, 1994. 377p.

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