Rota / Adeno IMC

Rota / Adeno IMC | Rota / Adeno IMC
Ref. 30.014.00

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Dr. Gilson Serio Pizzo

CRF MG - 5310

Anvisa 80027310317

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ENGLISH
BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to qualitatively determine Rotavirus and Adenovirus in human feces. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 30 °C.
- The product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed. Do not use products whose shelf life has expired.
- After opening the individual packaging, the device must be immediately used.

WORKING PRINCIPLE

Method: Immunochromatographic
The Rotavirus/Adenovirus rapid test is a lateral flow immunoassay for the qualitative detection of Rotavirus and Adenovirus antigens in human feces. During the test, the sample reacts with colloidal gold particles coated with anti-rotavirus and anti-adenovirus antibodies that migrates on the membrane by capillarity and react with a second anti-rotavirus antibody and a second anti-adenovirus antibody fixed in the test lines T1 and T2, respectively. The formation of these complexes generates a colored line in the T1 region indicating a positive result for rotavirus, and in the T2 region indicating a positive result for adenovirus. The absence of color in the test lines indicates a negative result. To validate the test, a colored line will appear in the control line region (C), indicating that an adequate sample volume was added and that the migration on the membrane occurred correctly.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: feces
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation: use recent collected feces samples. If necessary, store it at 2 to 8°C for, at maximum, 72 hours. For longer periods, store the samples at -20°C.

PRODUCT DESCRIPTION

DISP	One-time use dispositive with a nitrocellulose strand covered with anti-human Rotavirus and anti-human Adenovirus antibodies and goat anti-rat IgG antibodies.	
DIL	Phosphate buffer 5 g/L; NaCl 5 g/L; casein 3 g/L and preservative 0,02%.	

QUALITY CONTROL

Internal Control: the dispositive has an internal procedural control in the control line region (C). The appearance of a colored line in this region indicates that an adequate sample volume was added and that the sample migration occurred correctly.

External Control: it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure.

NECESSARY EQUIPMENT FOR TESTING

- Test tubes and clock or stopwatch.

TEST PROCEDURE AND INTERPRETATION

A) TEST PROCEDURE

Sample collection:

Collect approximately 1 to 2 mL or 1 to 2 g of feces in a clean and dry sample collection container. Preferably, perform the assay in up to 6 hours after sample collection.

Sample processing and test procedure:

1. Wait the reagents and samples to reach room temperature.
2. Process the sample as indicated below:

Solid Samples:


Unscrew the cap of the sample collection tube, then randomly stab the applicator into the fecal sample in at least 3 different sites to collect approximately 50 mg of feces. Insert the applicator into the tube with the extraction buffer, close the lid and mix it vigorously to homogenize the sample and the extraction buffer.

Liquid Samples:

With a dropper, aspirate the fecal sample and transfer two drops (approximately 50 µL) to the tube with the extraction buffer. Close the lid and mix it vigorously to homogenize the sample and the extraction buffer.

3. After homogenizing the sample with the extraction buffer, wait 2 minutes and apply two drops (approximately 80 µL) in the sample well (S) of the test dispositive. Start the chronometer. Avoid bubbles in the well sample.
4. After 10 minutes, read the results. Do not wait over 20 minutes to perform the reading.

If the sample does not migrate on the strand due to the presence of particles, centrifuge the sample homogenized with the extraction buffer, collect 80 µL of the supernatant and apply in the sample well of a new test dispositive, as indicated in the item 3 above.

B) ANALYSIS OF THE RESULTS

ROTAVIRUS POSITIVE: Two distinct colored lines appear,* one in the control line region (C) and another in the test line 1 region (T1).

ADENOVIRUS POSITIVE: Two distinct colored lines appear,* one in the control line region (C) and another in the test line 2 region (T2).

ROTAVIRUS AND ADENOVIRUS POSITIVE: Three distinct colored lines appear,* one in the control line region (C) and other two in the test line 1 and 2 regions (T1/T2).

* The intensity of the color in the test line regions (T1/T2) will vary depending on the concentration of rotavirus and adenovirus antigen present in the sample. Therefore, any shade of color in the test line region (T1/T2) should be considered positive.

NEGATIVE: One colored line appears in the control line region (C). No line appears in the test lines 1 and 2 regions (T1/T2).

INVALID: Control line (C) fails to appear. Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device.

C) INTERPRETATION

Rotavirus is the main agent responsible for acute gastroenteritis, especially in young children. Rotavirus is highly transmissible and is excreted in feces in high concentration. Transmission occurs through the fecal-oral route, from person to person or through contaminated water, food or objects, with an incubation period of 1 to 3 days. The clinical symptoms of the disease range from mild and short-lived diarrhea to severe cases of dehydration caused by profuse diarrhea and vomiting, accompanied by high fever. Rotaviral gastroenteritis can result in mortality for at-risk populations such as infants, elderly, and immunocompromised patients.

After rotavirus, enteric adenoviruses are the main cause of diarrhea in children. Although infection in children under two years old is more common, individuals of all ages can be affected. Adenoviruses usually cause respiratory illness, and the symptoms can range from the common cold to pneumonia, sore throat, and bronchitis. Depending on the type of adenovirus, other systems can be affected causing gastroenteritis, conjunctivitis, cystitis and, less commonly, neurological disease. Early diagnosis of diseases is extremely important to discriminate against other gastroenteritis and to guide medical conduct, minimizing the use of antibiotics.

PERFORMANCE CHARACTERISTICS

Clinical Sensitivity and Specificity:

Clinical sensitivity and specificity were evaluated by comparison with a reference kit of Latex Agglutination. It was obtained for Rotavirus and Adenovirus, respectively, a sensitivity of 97.3% and 95.2% and a specificity of 97.1% and 97.7%.

		Reference Kit		Total of Results
		Positive	Negative	
Rotavirus Kit Biotécnica	Positive	251	7	258
	Negative	7	236	243
Total of Results		258	243	501

		Reference Kit		Total of Results
		Positive	Negative	
Adenovirus Kit Biotécnica	Positive	118	6	124
	Negative	6	251	257
Total of Results		258	124	257

Analytical Sensitivity:

The limit of detection for Rotavirus and Adenovirus is of 1.0 x 10⁵ TCID₅₀/mL. Therefore, positive results indicate a viral concentration equal or bigger than 1.0 x 10⁵ TCID₅₀/mL.

Cross-Reactivity:

Cross reactivity was evaluated with the following organisms, not being detected cross-reaction up to the concentration of 1.0 x 10⁹ organisms/mL.

<i>Acinetobacter spp</i>	<i>Enterococcus faecium</i>	<i>Neisseria gonorrhea</i>
<i>Acinetobacter calcoaceticus</i>	<i>Gardnerella vaginalis</i>	<i>Enterococcus faecalis</i>
<i>Branhamella catarrhali</i>	<i>Streptococcus Grupo B</i>	<i>Proteus mirabilis</i>
<i>Salmonella choleraesius</i>	<i>Streptococcus Grupo C</i>	<i>Proteus vulgaris</i>
<i>Chlamydia trachomatis</i>	<i>Hemophilus influenzae</i>	<i>E.coli</i>
<i>Pseudomonas aeruginosa</i>	<i>Klebsiella pneumonia</i>	<i>Candida albicans</i>
<i>Neisseria meningitides</i>	<i>Staphylococcus aureus</i>	

Analytical Specificity:

Interfering substances up to the values presented below do not cause alterations in the results.

Bilirubin	Aspirin	Urea	Glucose	Caffeine
100 mg/dL	20 mg/dL	2000 mg/dL	2000 mg/dL	40 mg/dL

Ascorbic Acid	Oxalic Acid	Uric Acid	Albumin
20 mg/dL	60 mg/dL	60 mg/dL	2000 mg/dL

Precision:

Within-run precision was determined with 3 lots, one negative sample and 3 positive samples in 10 replicates for 10 days. Repeatability was observed in 100% of the replicates.
Between-run precision was determined with 3 lots, one negative sample and 3 positive samples in 10 replicates in one day. Reproducibility was observed in 100% of the replicates.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Do not replace components of this kit with those from other manufacturers, nor use components of different lots and codes.
- The test reading should be done in the period preconized in the test procedure to avoid false negatives.
- Concentrations below the detection limit may cause negative results.
- A negative result does not exclude the possibility of Rotavirus and/or Adenovirus infection. The result of this test should always be interpreted by a professional, considering the patient's clinical data. Failure to detect Rotavirus and/or Adenovirus by the test can be caused by several factors including not following some recommendations mentioned in this product insert, such as: storage condition, test procedure and sample collection, preparation and storage.
- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material. The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- Wilhelmi, I.; Roman, E.; Sanchez-Fauquier, A. Viruses causing gastroenteritis. Clin Microbiol Infect., v.9, p. 247-262, 2003.
- World Health Organization. Rotavirus: Vaccine Preventable Diseases Surveillance Standards (2018). Rotavirus. World Health Organization
- Cubitt, W. D. Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today., v. 1, p. 33-38, 1982.
- Nishio, Osamu, M. Ooseto, K. Takagi, Y. Yamasita, Y. Ishihara, and S. Isomura. Enzyme-Linked Immunosorbent Assay Employing Monoclonal Antibodies for Direct Identification of Enteric Adenoviruses (Ad40, 41) in Feces. Microbiol. Immunol. 1990; 34(10): 871-877.
- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens. Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
- CDC. National Center for Immunization and Respiratory Diseases, Division of Viral Diseases, 2019.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.

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