SANGUE OCULTO

FOB (Fecal Occult Blood) /
SANGRE OCULTA
Ref. 30.011.00

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ENGLISH

 BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the qualitative determination of human occult blood in stool. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 30 °C.
- The product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed.
- After opening the individual packaging the product is stable for 1h under the following humidity / temperature conditions: 35% / 20 ° C, 48% / 24,6 ° C and 85% / 37 ° C.
- Ready-to-use strip.

WORKING PRINCIPLE

Method: Immunochromatographic

Chromatographic immunoassay on nitrocellulose membrane. The human hemoglobin present in the sample binds to an antibody conjugated with colloidal gold and migrates through the membrane by capillarity to the test region containing human anti-hemoglobin antibodies forming a colored band. The unbound conjugate continues to migrate to the control region where a colored band should always appear.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: stool

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Stool	2 to 8 °C	3 days

PRODUCT DESCRIPTION

Device

Single-use nitrocellulose strip-containing casing coated with colloidal gold-conjugated monoclonal antibodies, anti-human hemoglobin monoclonal antibodies.

Diluent

Stabilizer.

QUALITY CONTROL

The device has a control line that must be visible for the interpretation of the result.

NECESSARY EQUIPMENT FOR TESTING

- Clock or stopwatch.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Before opening the packaging of the test device and the diluent, leave them at room temperature for 30 minutes (15 to 30 ° C).

Dilution of sample:

1. Insert the stick in at least 3 different locations in the sample.
2. Return the stick to the collection tube, close tightly and shake vigorously. Wait 2 minutes.

Performance of the test


1. Remove the device from its packaging and leave it on a horizontal surface.
2. Unscrew the cap from the diluent and add 2 drops (80 µL) of the diluted sample into the application scoop (S). Set the stopwatch.
3. Observe the result at exactly 5 minutes.

B) ANALYSIS OF RESULTS

Positive: two colored lines appear (red or pink), one line in the control region (C) and another in the test region (T).

Negative: one colored line appears only in the control region (C).

Invalid: no colored line appears in the control region (C).

C) INTERPRETATION

Investigation of bleeding in the gastrointestinal tract through stool occult blood research is an aid to clinical diagnosis and monitoring of bleeding disorders of the digestive tract. It is an indicated analysis for the preventive and anticipated screening of colorectal diseases in asymptomatic patients.

PERFORMANCE CHARACTERISTICS	
	Operating range
	Up to 4000 µg/mL
	Sensitivity - Detection Limit
	0.05 µg/mL
Analytical Specificity	
Sheep, chicken, ox, horse and pig hemoglobin	
500	

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	100
Positive Agreement	96%
Negative Agreement	100%

Prozone High Dose Effect:

There is no prozone effect up to 4000 µg/mL human hemoglobin.

Precision:

It was determined using a sample with concentration of 0.05 µg/mL of human hemoglobin in 20 replicates. The results obtained were in agreement with those expected.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Only use test and buffer devices whose packages are perfectly sealed.
- The test reading should not exceed 5 minutes.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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- MANDEL, Jack S. et al. Sensitivity, Specificity, and Positive Predictivity of the Hemoccult Test in Screening for Colorectal Cancers. Gastroenterology, v. 97, p.597-600, 1989.
- CATOMERIS, Peter et al. Effect of Temperature and Time on Fecal Hemoglobin Stability in 5 Fecal Immunochemical Test Methods and One Guaiac Method. Arch Pathol Lab Med, v. 142, p.75-82, jan. 2018.
- RABENECK, Linda et al. Fecal immunochemical tests compared with guaiac fecal occult blood tests for population-based colorectal cancer screenig. Can J Gastroenterol, v. 26, n. 3, p.131-147, mar. 2012.
- UCHIDA, Kazuo et al. Immunochemical detection of human blood in feces. Clinica Chimica Acta, v. 189, p.267-274, 1990.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.

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