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Enzymatic Sodium VET / Sodio Enzimatico VET  
Ref. 90.029.00

### FINALIDADE

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#### ENGLISH

#### INTENDED USE

Kit intended to determine sodium in serum samples. Veterinary diagnostic use only.


#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Enzymatic

Sodium is determined enzymatically through the activity of sodium-dependent β-galactosidase, which reacts with the ONPG (O-nitrophenyl-β-D-galactopyranoside) substrate to form O-nitrophenol and galactose. The formation rate of O-nitrophenol, which is spectrophotometrically measured at 405 nm, is proportional to the sodium concentration in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** Serum.

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum	4 to 8 °C	2 weeks
	-20 °C	1 year

#### PRODUCT DESCRIPTION

<b>R 1</b>	Tris buffer 450 mmol/L pH 9.0; Cryptand 5.4 mmol/L; β-Galactosidase ≥ 0.8 U/mL	<b>X</b>
<b>R 2</b>	Tris Buffer 10.0 mmol/L pH 9.0; α-nitrophenyl-β-D-galactopyranoside-ONPG 5.5 mmol/L.	<b>X</b>
<b>STD 1</b>	NaCl solution with sodium concentration indicated on the vial label.	<b>X</b>
<b>STD 2</b>	NaCl solution with sodium concentration indicated on the vial label.	<b>X</b>

Traceability of the standards STD1 and STD2 was obtained with the GBW09152 secondary standard of concentration defined by AAS.

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the controls below:

Normal Control – Quantinorm Vet	<b>REF</b>	90.040.00
Pathological Control – Quantialt Vet		90.041.00

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

1. Pipette into test tube:

	STD1	STD2	Sample
<b>R1</b>	720 µL	720 µL	720 µL
<b>Sample</b>	32 µL	32 µL	32 µL

2. Homogenize and incubate a 37 °C for 5 minutes.

3. Add:

<b>R2</b>	240 µL	240 µL	240 µL
-----------	--------	--------	--------

4.Homogenize and transfer immediately cuvette. Set the stopwatch. Enter the absorbance at 60 seconds (A<sub>1</sub>) and at 180 seconds (A<sub>2</sub>) at 405 nm by zeroing purified water.

#### B) Calculations

Sódio (mmol/L) = (Δ Sample Absorbance x Calibration Factor) – Intersection

Calibration Fator =  $\frac{[STD2]-[STD1]}{0,329 - 0,252}$

ΔAbs STD2 - ΔAbs STD1

Intersection = Calibration Fator x ΔAbs STD1 – [STD1]

Where:

Δ Sample Absorbance = (A<sub>2</sub> – A<sub>1</sub>) of the sample

[ STD2] = STD2 concentration

[ STD1] = STD1 concentration

Δ Abs STD2 = (A<sub>2</sub> – A<sub>1</sub>) of STD2

Δ Abs STD1 = (A<sub>2</sub> – A<sub>1</sub>) of STD1

**Example:**

		Absorbance		ΔAbs
		Abs 1	Abs 2	A2-A1
<b>STD1</b>	120 mmol/L	0,203	0,455	0,252
<b>STD2</b>	160 mmol/L	0,264	0,593	0,329
<b>Sample</b>		0,381	0,673	0,292

Calibration Factor =  $\frac{160,0 - 120,0}{0,329 - 0,252} = \frac{40}{0,077} = 519,5$

Intersection = (519,5 x 0,252) - 120,0 = 10,9

Sample concentration (mmol/L) = (0,292 x 519,5) - 10,9 = 140,8

**Automation:** This procedure can be applied on most automated analyzers. The protocols are available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br).

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The function of electrolytes in the human body is manifold. Among them are the maintenance of osmotic pressure and water distribution in the various liquid compartments of the organism, maintenance of adequate pH, regulation of the function of the heart and other muscles, participation in oxidation-reduction reactions and in catalysis as enzymatic cofactors. Abnormal electrolyte levels can be the cause or result of many disorders. Sodium is the main cation of the extracellular fluid. It plays a central role in maintaining normal water distribution and osmotic pressure. Decreased plasma sodium, hyponatremia, occurs as a consequence of decreased intake or excessive loss, this occurs in cases of intense sweating, prolonged vomiting, persistent diarrhea or enteropathies with salt loss. Kidney loss occurs from aldosterone and other mineralocorticoid deficiencies or from severe polyuria. Several other conditions are also the cause of hyponatremia: metabolic acidosis, edema, ascites due to chronic heart failure, uncontrolled diabetes, liver cirrhosis, nephrotic syndrome, malnutrition, inadequate secretion of antidiuretic hormone, etc. Hypernatremia occurs due to different clinical conditions, such as: hyperaldosteronism, increased production of mineralocorticoids, heart failure, liver disease, pregnancy, burns, osmotic diuresis, among others.

#### PERFORMANCE CHARACTERISTICS

Operating range
42.71 to 172.30 mmol/L

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
40.00 mmol/L

Analytical Specificity		
Species	Hemoglobin	Bilirubin
Canines	300 mg/dL	40 mg/dL
		2000 mg/dL

Felines	300 mg/dL	40 mg/dL	2000 mg/dL
Equinos	300 mg/dL	40 mg/dL	2000 mg/dL
Bovinos	500 mg/dL	40 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

	Accuracy			
	Canines	Felines	Equinos	Bovinos
Number of Samples	15 in duplicate	8 in duplicate	7 in duplicate	7 in duplicate
Regression Equation	$y = 1.10102x - 2.2606$	$y = 1.0896x - 13.61$	$y = 1.0897x - 11.981$	$y = 1.0346x - 5.1848$
Correlation Coefficient (R)	0.9757	0.9943	0.9815	0.9769

#### Precision:

Determined with two runs in duplicate.  
CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.

Samples (mmol/L)	Repetitions	Within-Run Precision		Samples (mmol/L)	Repetitions	Total Precision	
		SD (mg/dL)	CV (%)			SD (mg/dL)	CV (%)
68,39	20	2,7	4,0	68,08	40	2,16	3,2
156,98	20	2,22	1,4	157,31	40	2,37	1,5

- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the requirements by means of periodic checks.

#### REFERENCE RANGES

Canines	147 - 156 mmol/L
Felines	141 - 152 mmol/L
Equinos	132 - 152 mmol/L
Bovinos	132 - 146 mmol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units (mEq/L):

Sodium (mmol/L) x 1,0 = Sodium (mEq/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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