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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE


Kit intended to determine sodium in serum samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic
Sodium is determined enzymatically through the activity of sodium-dependent β-galactosidase, which reacts with the ONPG (O-nitrophenyl-β-D-galactopyranoside) substrate to form O-nitrophenol and galactose. The formation rate of O-nitrophenol, which is spectrophotometrically measured at 405 nm, is proportional to the sodium concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	2 weeks
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1	Tris buffer 450 mmol/L pH 9.0; Cryptand 5.4 mmol/L; β-Galactosidase ≥ 0.8 U/mL	X
R 2	Tris Buffer 10.0 mmol/L pH 9.0; α-nitrophenyl-β-D-galactopyranose - ONPG 5.5 mmol/L.	X
STD 1	NaCl solution with sodium concentration indicated on the vial label.	X
STD 2	NaCl solution with sodium concentration indicated on the vial label.	X

Traceability of the standards STD1 and STD2 was obtained with the GBW09152 secondary standard of concentration defined by AAS.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the controls below:

Normal Control – Quantinorm	REF	13.003.00
Pathological Control – Quantialt		13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Standard 1	Standard 2	Sample
STD	32 μL	32 μL	-
Sample	-	-	32 μL
R1	720 μL	720 μL	720 μL
Homogenize and incubate at 37 °C for 5 minutes			
R2	240 μL	240 μL	240 μL

2. Homogenize and incubate in a thermostated bucket. Trigger the timer.
3. Measure absorbance at 405 nm at 60 seconds (A1) and 180 seconds (A2), hitting zero with purified water.

B) CALCULATIONS

Sodium (mmol/L) = (Sample's ΔAbs x Calibration Factor) - Intersection

Calibration Factor = $\frac{\text{STD2 Concentration} - \text{STD1 Concentration (mmol/L)}}{\Delta\text{Abs STD2} - \Delta\text{Abs STD1}}$

Intersection = Calibration Factor x ΔAbs STD1 - STD1 Concentration

Where:

Sample Absorbance = (A2 – A1) of the sample

Absorbance STD2 = (A2 – A1) of STD2

Absorbance STD1 = (A2 – A1) of STD1

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaIda.ind.br

C) INTERPRETATION

The function of electrolytes in the human body is multiple. Among them are osmotic pressure and water distribution in the various fluid compartments of the organism, maintenance of adequate pH, regulation of heart and other muscle function, participation in oxidoreduction reactions and in catalysis as enzymatic cofactors. Abnormal levels of electrolytes may be the cause or the consequence of various disorders. Sodium is the major cation of the extracellular fluid. It plays a central role in maintaining the normal distribution of water and osmotic pressure. Decreased plasma sodium, hyponatremia, occurs as a consequence of decreased intake or excessive loss, which occurs in cases of severe sweating, prolonged vomiting, persistent diarrhea, or salt-losing enteropathies. Kidney failure occurs due to aldosterone deficiency and other mineralocorticoids or severe polyuria. Several other conditions also cause hyponatremia: metabolic acidosis, edema, ascites due to chronic heart failure, uncontrolled diabetes, liver cirrhosis, nephrotic syndrome, malnutrition, inadequate secretion of the antidiuretic hormone, etc. Hypernatremia occurs due to different clinical conditions, such as: hyperaldosteronism, increased mineralocorticoid production, heart failure, liver disease, pregnancy, burns, osmotic diuresis, among others.

PERFORMANCE CHARACTERISTICS

Operating range
63.0 to 180.0 mmol/L

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
5.80 mmol/L	63.33 mmol/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
500 mg/dL	40 mg/dL	1500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 0.993x - 0.030$
Correlation Coefficient (R)	0.9993

By applying the regression equation, the total systematic error estimated for decision levels of 70 mmol/L and 135 mmol/L was -0.65% and -0.68%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mmol/L)	Replicates	Within-Run Precision		Total Precision	
		SD (mmol/L)	CV (%)	SD (mmol/L)	CV (%)
82.90	80	0.68	0.80	0.71	0.90
142.90	80	0.99	0.70	1.27	0.90
174.50	80	1.46	0.80	1.77	1.00

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Serum	136 - 1145 mmol/L
-------	-------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Sodium (mmol/L) x 1.0 = Sodium (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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