# Transferrina

Transferrin | Transferrina  
Ref. 20.016.00

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### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

### INTENDED USE

Kit for the determination of transferrin in serum and plasma. *In vitro* diagnostic use.

### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- R1: Ready to use reagent.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

### WORKING PRINCIPLE

**Method:** Immunoturbidimetric anti-human transferrin antibodies form insoluble complexes when mixed with samples containing transferrin. The turbidity intensity is proportional to the transferrin concentration of the sample, which is spectrophotometrically measured at 540 nm.

### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum and plasma (EDTA and heparin)

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

### Preservation:

	Temperature	Stability Period
Serum and Plasma	2 to 8 °C	7 days
	-20 °C	3 months

### PRODUCT DESCRIPTION

**R 1**

Goat anti-human transferrin antibodies; activator, phosphate buffer pH 7.4; preserva-tive.

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Multicalibrador Plus	<b>REF</b>	21.011.00
Medium Multiparameter Control		21.003.00

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 540 nm.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnica.tda.ind.br](http://www.biotecnica.tda.ind.br)

REACTION PARAMETERS	
Type of reaction	Final Point / Growing
Wave length	540 nm
Temperature	37 °C
Sample volume/CAL*	3 µL
R1 volume*	300 µL
Incubation time	4 minutes
Calibration	Calibration curve**

\* Sample/CAL and R1 volumes can be modified as long as the established ratio is maintained.

### \*\*CALIBRATION CURVE:

Dilute the PLUS Multicalibrador in 9g / L NaCl as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	10	25	50	75	100
NaCl 9g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1

Multiply the concentration of transferrin indicated in the instruction of use of the Multicalibrador PLUS by the corresponding factor to obtain the transferrin concentration of each dilution. The transferrin concentration of the calibrator is traceable to Standard Reference Materials ERM DA470k / IFCC.

### B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

### C) INTERPRETATION

Transferrin is the major plasma protein for transporting iron and is synthesized in the liver. The transferrin molecule binds specifically to Fe<sup>3+</sup>, forming the transferrin-Fe<sup>3+</sup> complex, which circulates freely through the serum transport-ing the iron to the storage sites. The evaluation of transferrin levels is useful for differential diagnosis of anemia, to monitor its treatment and to assess the nutritional status of the patient. The level of transferrin increases in hypo-chromic anemia (iron deficiency), in pregnancy and during treatment with estrogens, and decreases in inflammation and malignant tumors.

### PERFORMANCE CHARACTERISTICS

Operating range
15,774 to 606 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mL (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.031 mg/dL	15.774 mg/dL

Analytical Specificity			
Hemoglobin	Bilirubin	Triglycerides	Rheumatoid Factor
500 mg/dL	40 mg/dL	1000 mg/dL	170 IU/mL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	53 in duplicate
Regression Equation	y = 0.982x - 2.2593
Correlation Coefficient (R)	0.9970

By applying the regression equation, the total systematic error estimated for decision levels of 200 mg/dL and 360 mg/dL was -2.90% and -2.40%, respectively.

### Prozone High Dose Effect:

There is no prozone effect up to 1541 mg/dL transferrin.

### Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
50.67	80	0.931	1.8	1.352	2.7
359.73	80	5.534	1.5	10.548	2.9
441.64	80	7.837	1.8	9.983	2.3

CV: Coefficient of variation; SD: Standard deviation

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

### REFERENCE RANGES

Newborn	117 to 250 mg/dL
Adults (20-60 years)	200 to 360 mg/dL
Adults (> 60 years)	160 to 340 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling +55 35 3214 4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).

To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).

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