Triglicérides VET

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ENGLISH

INTENDED USE

Kit intended to determine triglycerides in serum and plasma samples. Veterinary diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic - Colorimetric

Triglycerides are hydrolyzed by lipoprotein lipase to form free glycerol, which is transformed into glycerol-3-phosphate by the action of glycerol kinase. The glycerol 3-phosphate is oxidized to dihydroxyacetone and hydrogen peroxide by glycerol-P-oxidase. Hydrogen peroxide reacts with 4-aminoantipyrine and 4-chlorophenol under the catalytic action of peroxidase, producing a pink compound (quinonimine) that can be spectrophotometrically measured at 505 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (EDTA).

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	7 days
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1	Pipes buffer ≥ 20 mmol/L; 4-chlorophenol ≥ 1 mmol/L; 4-aminoantipyrine ≥ 0,1 mmol/L; adenosine triphosphate ≥ 0,5 mmol/L; glycerol kinase ≥ 500 U/L; peroxidase ≥ 1000 U/L; lipoprotein lipase ≥ 1000 U/L; glycerol-3-phosphate oxidase ≥ 1000 U/L; activators; detergents, stabilizers and preservative.	
	Glycerol is equivalent to the triglyceride concentration of 200 mg/dL, preservatives. Traceable to reference material NIST 1951b.	

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal Vet		90.039.00
Normal Control – Quantinorm Vet	REF	90.040.00
Pathological Control – Quantialt Vet		90.041.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 505 nm (490 – 510 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 10 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 505 nm (490 – 510 nm). The final reaction is stable for 25 minutes.

B) CALCULATIONS

Triglycerides (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Triglycerides (mg/dL) = Sample Absorbance x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaitda.ind.br

C) INTERPRETATION

Triglycerides are esters of glycerol. Part is synthesized in the liver and another part is obtained from the diet. Through the action of pancreatic and intestinal lipases in the presence of bile acids, they are first hydrolyzed to glycerol, monoglycerides and fatty acids. After absorption, these components again form triglycerides in the intestinal epithelial cells and then are combined with cholesterol and apolipoproteins to form chylomicrons, which are secreted into lymphatic system to then reach the circulation. Triglycerides are the primary metabolic fuel carried by chylomicrons and are distributed to the liver and peripheral cells, where they are hydrolyzed into fatty acids by lipases. The determination of triglycerides is used in the diagnosis and treatment of animals with diabetes mellitus, nephrosis, liver obstruction, disorders in the metabolism of lipids and several other endocrine diseases.

PERFORMANCE CHARACTERISTICS

Operating range
9.34 to 801.93 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
9.00 mg/dL

Analytical Specificity		
Species	Hemoglobin	Bilirubin
Caninos	300 mg/dL	40 mg/dL
Felinos	200 mg/dL	10 mg/dL
Equinos	200 mg/dL	10 mg/dL
Bovinos	100 mg/dL	10 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy – Canines	
Number of Samples	40 in duplicate
Regression Equation	$y = 0,9775x + 0,8771$
Correlation Coefficient (R)	0,9974
Accuracy – Felines	
Number of Samples	20 in duplicate
Regression Equation	$y = 1,0077x - 2,4532$
Correlation Coefficient (R)	0,9989
Accuracy – Equines	
Number of Samples	11 in duplicate
Regression Equation	$y = 0,98x - 1,0945$
Correlation Coefficient (R)	0,9979
Accuracy – Bovines	
Number of Samples	10 in duplicate
Regression Equation	$y = 0,961x + 0,0824$
Correlation Coefficient (R)	0,9951

Precision:

Determined with two runs in duplicate.

Samples (mg/dL)	Repetitions	Within-Run Precision		Samples (mg/dL)	Repetitions	Total Precision	
		SD (mg/dL)	CV (%)			SD (mg/dL)	CV (%)
13,21	20	0,61	4,6	12,99	40	0,65	5,0
190,28	20	5,31	2,8	183,93	40	8,88	4,8

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Canine	20 – 112 mg/dL
Feline	10 – 114 mg/dL
Equine	4 – 44 mg/dL
Bovine	0 – 14 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Triglycerides (mg/dL) x 0,0113 = Triglycerides (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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