Analytical Specificity	
Hemoglobin	Bilirubin
500 mg/dL	35 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 0.995x - 0.950$
Correlation Coefficient (R)	0.9997

By applying the regression equation, the total systematic error estimated for decision levels of 180 mg/dL and 500 mg/dL was -1.02% and -0.69%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
104.69	80	1.173	1.10	1.176	1.10
227.77	80	3.212	1.40	3.770	1.70
626.55	80	2.090	0.30	2.228	0.40

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Age	Triglyceride Values (mg/dL)	Interpretation
2 to 19 years	< 100	Desirable
	100 – 129	Limitrophe
	≥ 130	High
≥ 20 years	< 150	Desirable
	150 – 200	Limitrophe
	200 – 499	High
	≥ 500	Very High

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Triglycerides (mg/dL) x 0,0113 = Triglycerides (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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