Troponina I IMC

Troponin I IMC / Troponina I IMC
 Ref. 30.026.00

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Dr. Gilson Serio Pizzo
 CRF MG - 5310

Anvisa 80027310323

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to qualitatively determine Troponin I in serum, plasma or whole blood samples. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 30 °C.
- The product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed. Do not use products whose shelf life has expired.
- After opening the individual packaging, the device must be immediately used.
- Do not use the diluent after 6 months of opening the bottle.

WORKING PRINCIPLE

Method: Immunochromatographic

The Troponin I Rapid Test is a qualitative lateral flow immunoassay for the detection of Troponin I in whole blood, serum or plasma samples. During the test, the sample with Troponin I reacts with particles coated with an anti-Troponin I antibody, forming a complex that migrates through the membrane by capillarity. This complex reacts with a second anti-troponin I antibody attached to the test line, generating a colored line in this region that indicates a positive result. The absence of staining in the test line region indicates a negative result. To validate the test, a colored line will appear in the control line region, indicating that an adequate volume of sample has been added and that migration through the membrane has occurred correctly.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum, plasma (EDTA-K2, Sodium Heparin, Sodium Citrate and Sodium Oxalate)

and whole blood

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

| | |
|------------------|---|
| Serum and Plasma | Store it at 2 to 8°C for, at maximum, 72 hours. For longer periods, store the samples at -20°C. |
| Venous blood | Store it at 2 to 8°C for, at maximum, 24 hours. Do not freeze. |
| Capillary blood | Must be used immediately. Do not store or freeze. |

PRODUCT DESCRIPTION

One-time use dispositive with a nitrocellulose strand covered with anti-Troponin I antibodies, capture antibodies and goat anti-rat IgG antibodies.

DISP

DIL

Phosphate buffer 5 g/L; NaCl 5 g/L; casein 3 g/L and preservative 0,02%.

QUALITY CONTROL

Internal Control: the dispositive has an internal procedural control in the control line region (C). The appearance of a colored line in this region indicates that an adequate sample volume was added and that the sample migration occurred correctly.

External Control: it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure.

NECESSARY EQUIPMENT FOR TESTING

- Test tubes and clock or stopwatch.

TEST PROCEDURE AND INTERPRETATION

A) TEST PROCEDURE

Collection and sample processing:

- Capillary blood:

Wash the patient's hand with soap and warm water or clean the area with alcohol and let it dry. Rub down the hand towards the fingertip of the

middle or ring finger. Puncture the skin with a sterile lancet and wipe away the first sign of blood. Gently rub the hand from wrist to finger to form a round drop of blood over the puncture site.

- Serum and Plasma:

Serum and Plasma samples should be separated as soon as possible to avoid hemolysis.

Test procedure:

1. Wait the reagents and samples to reach room temperature.
2. Apply the sample as indicated below, avoiding bubble formation in the sample well.

- Serum or Plasma samples:

With a dropper, apply 3 drops of the sample (approximately 75 µL) in the sample well (S) of the test dispositive.


- Venous blood samples:

With a dropper, apply 3 drops of the sample (approximately 75 µL) in the sample well (S) of the test dispositive and then, apply 1 drop of buffer (approximately 40 µL).

- Capillary blood samples:

Fill the capillary tube with the blood sample from the fingertip (approximately 75 µL), transfer all the sample volume to the sample well of the dispositive (S) and then, apply 1 drop of buffer (approximately 40 µL).

3. Start the chronometer. After 10 minutes, read the results. Do not wait over 20 minutes to perform the reading.

B) ANALYSIS OF THE RESULTS

POSITIVE: Two distinct colored lines appear,* one in the control line region (C) and another in the test line region (T).

* The intensity of the color in the test line region (T) will vary depending on the concentration of Troponin I present in the sample. Therefore, any shade of color in the test line region (T) should be considered positive.

NEGATIVE: One colored line appears in the control line region (C). No line appears in the test line region (T).

INVALID: Control line (C) fails to appear. Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control

line failure. Review the procedure and repeat the test with a new test dispositive.

C) INTERPRETATION

The troponins are a complex of three proteins (troponins C, T, and I) that regulate muscle contraction in cardiac and skeletal muscle through the calcium-dependent interaction of myosin with actin. Troponin I has an isoform expressed only in cardiac muscle, which makes it a highly specific marker for myocardial injury. After myocardial injury, troponin I is released into the bloodstream, so its concentration can be detected within 3 to 4 hours. The pattern of troponin I release into the blood stream is similar to that of CK-MB, but the concentration of CK-MB returns to normal after 72 hours of myocardial injury, while troponin I values remain detectable for 6-10 days after myocardial injury, providing a larger window of detection for the diagnosis of myocardial injury. Troponin I can be used both as a marker of acute myocardial infarction and for the differential diagnosis of other pathologies that affect the heart tissue, such as: congestive heart failure, acute coronary syndrome, myocarditis and perioperative cardiac damage.

PERFORMANCE CHARACTERISTICS

Clinical Sensitivity and Specificity:

Clinical sensitivity and specificity were evaluated by comparison of 445 samples with a reference kit. It was obtained a sensitivity of 97.6% and a specificity of 99.4%.

| | | Reference (Chemiluminescence) | | Total of Results |
|------------------|----------|-------------------------------|----------|------------------|
| | | Positive | Negative | |
| Biotécnica | Positive | 83 | 2 | 85 |
| | Negative | 2 | 358 | 360 |
| Total of Results | | 85 | 360 | 445 |

Analytical Sensitivity:

The limit of detection is of 0.5 ng/mL. Therefore, positive results indicate a Troponin I concentration equal or bigger than 0.5 ng/mL.

Cross-Reactivity:

Cross reactivity was evaluated with positive samples for other muscular protein and diseases, not being detected cross-reaction with:

| | | |
|---------------------|-----------------|-----------------------|
| Hepatitis B (HBcAb) | HIV | Rheumatoid factor |
| Hepatitis B (HBsAg) | <i>H.pylori</i> | Mononucleosis (IgM) |
| Hepatitis B (HBsAb) | Syphilis | Toxoplasmosis (IgG) |
| Hepatitis B (HBeAg) | Rubella (IgG) | Cytomegalovirus (IgG) |
| Hepatitis B HBeAb) | Troponin T | Cardiac myosin |
| Skeletal Troponin I | | |

Analytical Specificity:

Interfering substances up to the values presented below do not cause alterations in the results.

| Ascorbic Acid | Gentisic Acid | Oxalic Acid | Acetylsalicylic Acid |
|---------------|---------------|-------------|----------------------|
| 20 mg/dL | 20 mg/dL | 600 mg/dL | 20 mg/dL |

| Bilirubin | Hemoglobin | Creatin | Acetaminophen |
|------------|------------|-----------|---------------|
| 1000 mg/dL | 1000 mg/dL | 200 mg/dL | 20 mg/dL |

| Cholesterol | Albumin | Caffeine | Triglycerides |
|-------------|-------------|----------|---------------|
| 800 mg/dL | 10500 mg/dL | 20 mg/dL | 1600 mg/dL |

Precision:

Within-run precision was determined with one lot, one negative sample and 6 positive samples in 3 replicates for 10 days. Repeatability was observed in 100% of the replicates.

Between-run precision was determined with 3 lots, one negative sample and 6 positive samples in 3 replicates in 3 days. Reproducibility was observed in 100% of the replicates.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Do not replace components of this kit with those from other manufacturers, nor use components of different lots and codes.
- The test reading should be done in the period preconized in the test procedure to avoid false negatives.
- Concentrations below the detection limit may cause negative results.
- The result of this test should always be interpreted by a professional, considering the patient's clinical data. Failure to detect Troponin I by the test can be caused by several factors including not following some recommendations mentioned in this product insert, such as: storage condition, test procedure and sample collection, preparation, and storage.
- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material. The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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
- APPLE, Fred S. et al. Cardiac troponin assays: guide to understanding analytical characteristics and their impact on clinical care. Clinical chemistry, v. 63, n. 1, p. 73-81, 2017.
- JAVED, Usman et al. Frequency of elevated troponin I and diagnosis of acute myocardial infarction. The American journal of cardiology, v. 104, n. 1, p. 9-13, 2009.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.

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