Ureia UV VET

Urea UV VET / Urea UV VET
Ref. 90.011.00

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| Espécies | Intervalo de Referência |
|----------|-------------------------|
| Caninos | 10,0 – 28,0 mg/dL |
| Felinos | 20,0 – 30,0 mg/dL |
| Equinos | 10,0 – 45,0 mg/dL |
| Bovinos | 20,0 – 30,0 mg/dL |

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ENGLISH

INTENDED USE

Kit intended to determine urea in serum and plasma samples. Veterinary diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.
- **Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 4 weeks, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C). An absorbance lower than 0,900 for the reagent, in a spectrophotometer set to zero with purified water at 340 nm, indicates its deterioration.

WORKING PRINCIPLE

Method: Enzymatic UV

The urea from sample is hydrolyzed by the enzyme urease producing carbon dioxide and ammonium ions. The ammonium ions, in presence of the enzyme glutamate dehydrogenase and the substrates NADH and α -ketoglutarate, produces NAD⁺ and glutamate. The consumption rate of NADH in the assay can be measured spectrophotometrically measured at 340 nm, being proportional to the urea concentration.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum and plasma.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

| | Temperature | Stability Period |
|-------------------|-------------|------------------|
| Serum and Plasma* | 4 to 8 °C | 5 days |
| | -20 °C | 6 months |

*The anticoagulants citrate, fluoride and sodium oxalate interfere in the dosage of urea.

PRODUCT DESCRIPTION

| | | |
|------------|---|--|
| R 1 | TRIS buffe ≥ 100 mmol/L; EDTA ≥ 1 mmol/L; α -ketoglutarate ≥ 2 mmol/L; ADPK ≥ 1 mmol/L; urease ≥ 5000 U/L; glutamate dehydrogenase ≥ 1000 U/L; activators; detergent; stabilizers; preservative. | |
| R 2 | Carbonate buffer ≥ 2 mmol/L; NADH $\geq 0,5$ mmol/L; preservatives. | |
| STD | Phosphate buffer ≥ 2 mmol/L; preservative; urea in a concentration equivalent to 70 mg/dL. Traceable to reference material NIST 912a. | |

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

| | |
|--------------------------------------|----------------------|
| Calibrator - Autocal Vet | 90.039.00 |
| Normal Control – Quantinorm Vet | REF 90.040.00 |
| Pathological Control – Quantialt Vet | 90.041.00 |

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Heat the work reagent at 37 °C for 3 minutes.
2. Zero the reading equipment with purified water at 340 nm.
3. Pipette in the test tubes:

| | STD | Sample |
|--------------|------------|-------------|
| STD | 10 μ L | - |
| Sample | - | 100 μ L |
| Work Reagent | 1,0 mL | 1,0 mL |

4. Homogenize and insert it immediately in the thermostatic cuvette at 37°C. Activate the chronometer.

5. Measure the absorbance at 340 nm at 30 seconds (A1) and at 120 seconds (A2) of sample and standard.

B) CALCULATIONS

Urea UV (mg/dL) = $\frac{\text{Sample's } (A1 - A2)}{\text{STD'} (A1 - A2)} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD'} (A1 - A2)}$

Urea UV (mg/dL) = Sample's (A1 - A2) x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaltda.ind.br

C) INTERPRETATION

Urea is synthesized in the liver through arginase and is the main product of protein catabolism. Urea is filtered through the glomerular filtration system in a concentration equal to that of blood. When there is a higher flow rate, there is less urea absorption and vice versa. In situations where there is low glomerular filtration, the highest concentration of urea is observed. This causes an increase in blood concentration. Urea concentration is altered by extrarenal factors such as increased and prolonged protein intake. Due to these interferences, urea is not a good marker of renal function on an individual basis. To analyze renal function, this parameter must be adjusted to urinary creatine, protein and density levels. The reduction in urea levels may occur as a result of reduced production, such as in cases of liver failure, in cirrhosis, without porto-sinus shunting, and in cases of reduced dietary protein and hypoproteinemia.

PERFORMANCE CHARACTERISTICS

| Operating range |
|-----------------------|
| 10.17 to 378.83 mg/dL |

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

| Sensitivity |
|-------------|
| 10,00 mg/dL |

| Quantification Limit |
|----------------------|
| 10.00 mg/dL |

| Analytical Specificity | | | |
|------------------------|------------|-----------|---------------|
| Species | Hemoglobin | Bilirubin | Triglycerides |
| Canines | 500 mg/dL | 40 mg/dL | 2000 mg/dL |
| Felinos | 500 mg/dL | 40 mg/dL | 2000 mg/dL |
| Equinos | 500 mg/dL | 10 mg/dL | 2000 mg/dL |
| Bovines | 500 mg/dL | 40 mg/dL | 2000 mg/dL |

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy - Canines | |
|-----------------------------|---------------------|
| Number of Samples | 40 in duplicate |
| Regression Equation | y = 0.996x - 0.202 |
| Correlation Coefficient (R) | 0.9986 |
| Accuracy - Felinos | |
| Number of Samples | 20 in duplicate |
| Regression Equation | y = 1.002x - 0.549 |
| Correlation Coefficient (R) | 0.9942 |
| Accuracy - Equines | |
| Number of Samples | 11 in duplicate |
| Regression Equation | y = 0.963x + 1.215 |
| Correlation Coefficient (R) | 0.9633 |
| Accuracy - Bovines | |
| Number of Samples | 10 in duplicate |
| Regression Equation | y = 1.065x - 0.5328 |
| Correlation Coefficient (R) | 0.9993 |

Precision:

Determined with two runs in 20 replicates

| Samples (U/L) | Number of samples | Within-Run Precision | | Samples (U/L) | Number of samples | Total Precision | |
|---------------|-------------------|----------------------|--------|---------------|-------------------|-----------------|--------|
| | | SD (U/L) | CV (%) | | | SD (U/L) | CV (%) |
| 72.85 | 20 | 2.37 | 3.30 | 71.75 | 40 | 2.06 | 2.90 |
| 183.89 | 20 | 5.20 | 2.80 | 186.56 | 40 | 9.02 | 4.80 |

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

| Species | 37°C |
|---------|-----------------|
| Canine | 10.0 – 28.0 U/L |
| Feline | 20.0 – 30.0 U/L |
| Equine | 10.0 – 45.0 U/L |
| Bovine | 20.0 – 30.0 U/L |

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Urea (mg/dL) x 0,167 = Urea (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnicas products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnicas instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.com.br.

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- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACCPress, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p. 493-501, 1981.
- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. **Clinical biochemistry of domestic animals**. 5th. ed. San Diego: Academic Press, c1997 932p.
- THRALL, Mary Anna. **Hematologia e bioquímica clínica veterinária**. São Paulo: Roca, 2007. 582 p.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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