## Ureia Enzimática

Enzymatic Urea / Ureia Enzimática  
Ref. 10.013.00

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#### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

#### INTENDED USE

Kit intended to determine urea in serum, plasma and urine samples. Diagnostic use only.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.
- Work Reagent:** mix in the proportion of 25 parts of R1 + 1 part of R3. Homogenize gently. The work reagent is stable for 21 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).

#### WORKING PRINCIPLE

**Method:** Colorimetric enzymatic

The urea of the sample is hydrolyzed by urease with production of carbonic gas and ammonium ions. It reacts with salicylate and sodium hypochlorite in alkaline medium, in the presence of nitroprusside, producing indophenol, green in color, which can be measured in a spectrophotometer at 580 nm. The intensity of the color is proportional to the concentration of urea in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

**Sample Type:** serum, plasma (EDTA and heparin) and urine.

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

Temperature	Stability Period
-------------	------------------

Serum	4 to 8 °C	7 days
	-20 °C	1 year
24 Hours Urine*	4 to 8 °C	7 days
	-20 °C	28 days

#### Preparation:

**Urine:** work with samples collected in the period of 24 hours in a flask with 2 mL of HCI 50% v/v. Centrifuge at 3000 rpm for 10 minutes and collect the supernatant. Perform a 1:50 dilution of the supernatant with purified water and perform the assay. Multiply the result obtained by 50. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Phosphate buffer ≥ 3 mmol/L; Sodium salicylate ≥ 10 mmol/L; Sodium nitroprusside ≥ 1 mmol/L; stabilizers; preservative.	
<b>R 2</b>	Sodium hypochlorite ≥ 0,01% v/v; Sodium hydroxide ≥ 100 mmol/L.	
<b>R 3</b>	Phosphate buffer ≥ 10 mmol/L; Glycerol ≥ 20% v/v; Urease ≥ 50.000 U/L; stabilizers; preservative.	 
<b>STD</b>	Phosphate buffer ≥ 2 mmol/L; Urea at a concentration equivalent to 70 mg/dL; preservative. Traceable to NIST reference material 912a.	

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00
Normal Control - Quantinorm	13.003.00
Pathological Control - Quantial	13.004.00

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 580 nm (570 a 650 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	STD	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
Purified Water	10 µL	-	-
Working Reagent	1,0 mL	1,0 mL	1,0 mL
Homogenize and incubate for 5 minutes at 37 °C			
R2	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 5 minutes.

3. Measure Standard Absorbance (STD) and Sample vs. Blank at 580 nm. Color is stable for 15 minutes.

##### B) CALCULATIONS

$$\text{Ureia (mg / dL)} = \frac{\text{Sample Absorbance} \times \text{STD Concentration (mg / dL)}}{\text{Absorbance of STD}}$$

**Exemplo:**

STD concentration = 70 mg / dL

Sample Absorbance = 0.209

Absorbance of STD = 0.378

$$\text{Ureia (mg / dL)} = \frac{0,209}{0,378} \times 70 = 38,7 \text{ mg / dL}$$

**Calculations with the Calibration Factor:**

Calibration Factor =  $\frac{\text{STD Concentration (mg/dL)}}{\text{Absorbance of STD}}$

Ureia (mg/dL) = Sample Absorbance x Calibration Factor

**Urine:**

Determine the concentration of urea in the urine sample using the same calculation procedure for serum. Multiply the value obtained by 50 to obtain the result in mg / dL. This should be used to calculate the value of urea in g / 24h, according to the equation below:

$$\text{Ureia (g/24 hours)} = \frac{\text{Ureia (mg/dL)} \times \text{Urinary Volume (mL)}}{100}$$

\*The supernatant from the urine sample must be diluted in a 1:50 proportion. To determine the urea concentration in the urine sample (mg/dL), multiply the obtained result by 50.

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicalda.ind.br](http://www.biotechnicalda.ind.br)

#### C) INTERPRETATION

Urea is the main nitrogenous metabolite of protein catabolism. More than 90% of urea is excreted through the kidneys; therefore, renal disease is associated with accumulation of urea in the blood. The dosage of blood and urinary urea can provide useful clinical information. It can be elevated due to various factors, such as high protein diet, high protein catabolism, resorption of blood proteins after gastrointestinal bleeding, treatment with cortisol and its analogs, dehydration and reduced perfusion of the kidneys. It is also

elevated in obstructive post-renal conditions such as malignant tumors, nephrolithiasis and prostatism.

#### PERFORMANCE CHARACTERISTICS

	Operating range
Serum	2.206 to 200 mg/dL
Urine*	2.303 to 100 g/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

\*Corrected by dilution factor.

	Sensitivity	
	Detection Limit	Quantification Limit
Serum	0.770 mg/dL	2.206 mg/dL
Urine*	0.006 g/L	2.303 g/L

\*Corrected by dilution factor.

	Analytical Specificity		
	Hemoglobin	Bilirubin	Triglycerides
Serum	500 mg/dL	12 mg/dL	700 mg/dL
Urine	150 mg/dL	20 mg/dL	-

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine
Number of Samples	40 in duplicate	20 in duplicate
Regression Equation	y = 0.995x + 0.08	y = 1.013x + 0.231
Correlation Coefficient (R)	0.9997	0.9991

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 25 mg/dL and 150 mg/dL was -0.18% and -0.45%, respectively.

**Urine:** by applying the regression equation, the total systematic error estimated for decision levels 30 g/24h and 100 g/24h was 2,07% and 1,53%, respectively.

#### Precision:

Determined with two runs in duplicate per day for 20 days.

Serum Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
37,02	80	0,53	1,5	1,06	2,9
88,00	80	0,26	0,3	1,42	1,6
162,68	80	0,81	0,5	0,96	0,6

CV: Coefficient of variation; SD: Standard deviation

Urine Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
953,75	80	10,98	1,2	12,37	1,3
2676,03	80	20,79	0,8	21,98	0,8

CV: Coefficient of variation; SD: Standard deviation

#### RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

	Serum and Plasma	13 - 45 mg/dL	2.17 - 7.51 mmol/L
Urine		26 - 43 g/24 hours	

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Urea (mg/dL) x 0,167 = Urea (mmol/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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