	
Ureia UV	Responsável Técnico: Dr. Gilson Serio Pizzo CRF MG – 5310
Urea UV / Urea UV Ref. 10.012.00 / FD 38.053.00 / SD 39.062.00	Anvisa 80027310263

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ENGLISH

 **BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

INTENDED USE

Kit intended to determine urea in serum, plasma and urine samples. Diagnostic use only.



STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.
- **Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 4 weeks, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C). An absorbance lower than 0,900 for the reagent, in a spectrophotometer set to zero with purified water at 340 nm, indicates its deterioration.
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

WORKING PRINCIPLE

The urea from sample is hydrolyzed by the enzyme urease producing carbon dioxide and ammonium ions. The ammonium ions, in presence of the enzyme glutamate

dehydrogenase and the substrates NADH and α-ketoglutarate, produces NAD⁺ and glutamate. The consumption rate of NADH in the assay can be measured spectrophotometrically measured at 340 nm, being proportional to the urea concentration.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum, plasma and urine.
Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

	Temperature	Sample Stability
Serum and plasma*	4 to 8 °C	7 days
	-20 °C	1 year
24 Hours Urine	4 to 8 °C	7 days
	-20 °C	28 days

*The anticoagulants citrate, fluoride and sodium oxalate interfere in the dosage of urea.

Preparation:
Urine: work with samples collected in the period of 24 hours in a flask with 2 mL of HCl 50%. Centrifuge at 3000 rpm for 10 minutes and collect the supernatant. Perform a 1:50 dilution of the supernatant with purified water and perform the assay. Multiply the result obtained by 50. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

PRODUCT DESCRIPTION

R 1	TRIS buffe ≥ 100 mmol/L; EDTA ≥ 1 mmol/L; α-ketoglutarate ≥ 2 mmol/L; ADPK ≥ 1 mmol/L; urease ≥ 5000 U/L; glutamate dehydrogenase ≥ 1000 U/L; activators; detergent; stabilizers; preservative.	 
R 2	Carbonate buffer ≥ 2 mmol/L; NADH ≥ 0,5 mmol/L; preservatives.	
STD	Phosphate buffer ≥ 2 mmol/L; preservative; urea in a concentration equivalent to 70 mg/dL. Traceable to reference material NIST 912a.	

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Autocal H	13.002.00	
Normal Control – Quantinorm	13.003.00	
Pathological Control – Quantialt	13.004.00	

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

TEST PROCEDURE

1. Heat the work reagent at 37 °C for 3 minutes.
2. Zero the reading equipment with purified water at 340 nm.
3. Pipette in the test tubes:

	STD	Sample
STD	10 µL	-
Sample	-	10 µL
Work Reagent	1.0 mL	1.0 mL

3. Homogenize and insert it immediately in the thermostated cuvette at 37 °C. Activate the chromometer.
4. Measure the absorbance at 340 nm at 30 seconds (A1) and at 120 seconds (A2) of sample and standard.

B) CALCULATIONS

Calculation of Serum and Plasma
Urea UV (mg/dL) = $\frac{\text{Sample's (A1 - A2)}}{\text{STD' (A1 - A2)}}$ x STD Concentration (mg/dL)

Calculations with the Calibration Factor:
Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD' (A1 - A2)}}$

Urea UV (mg/dL) = Sample's (A1 - A2), x Calibration Factor
Calculations for urine samples:
Urea UV (mg/24 hours) = Urea UV (mg/dL) * x Urinary Volume (mL)
100

*The supernatant from the urine sample must be diluted in a 1:50 proportion. To determine the urea concentration in the urine sample (mg/dL), multiply the obtained result by 50.

Automaton: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Urea is the main nitrogenous metabolite of protein catabolism. More than 90% of urea is excreted through the kidneys; therefore, renal disease is associated with accumulation of urea in the blood. The dosage of blood and urinary urea can provide useful clinical information. It can be elevated due to various factors, such as high protein diet, high protein catabolism, resorption of blood proteins after gastrointestinal bleeding, treatment with cortisol and its analogs, dehydration and reduced perfusion of the kidneys. It is also elevated in obstructive post-renal conditions such as malignant tumors, nephrolithiasis and prostatism.

PERFORMANCE CHARACTERISTICS

Operating range	
13.3 to 200.0 mg/dL	

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
3.52 mg/dL	13.30 mg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
500 mg/dL	30 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 1.016x - 0.572
Correlation Coefficient (R)	0.9971

By applying the regression equation, the total systematic error estimated for decision levels of 50 mg/dL and 130 mg/dL was 0.46% and 1.16%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
31.73	80	0.346	1.10	0.474	1.50
108.31	80	0.801	0.70	1.654	1.50
198.70	80	1.345	0.70	2.764	1.40

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

Serum and Plasma	13 - 45 mg/dL	2.17 - 7.51 mmol/L
Urine	26 - 43 g/24 hours	

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Urea (mg/dL) x 0,167 = Urea (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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