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**ENGLISH**

 **BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

**INTENDED USE**  
Kit for the determination Rheumatoid Factor (RF) in serum. *In vitro* diagnostic use.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable for 2 months, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagents and controls.

**WORKING PRINCIPLE**  
**Method:** Slide haemagglutination  
Waaler Rose is a plate haemagglutination method for qualitative and semi-quantitative detection of rheumatoid factor. When the latter is present in the sample, it reacts with sheep red cells coated with purified rabbit IgG, which are agglutinated.

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**  
**Sample Type:** serum  
**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.  
**Preservation:**

	Temperature	Stability Period
Serum	20 to 25°C	1 day
	4 to 8 °C	8 days
	-20 °C	3 months

**PRODUCT DESCRIPTION**

**R 1** Suspension of erythrocytes sensitized with rabbit serum IgG anti-sheep erythrocytes pH 8,2; preservative.

**CONTROL +** Human serum with rheumatoid factor concentration > 30 IU/mL; preservative. (INCLUDED IN THE Waaler Rose KIT)

**CONTROL -** Human serum in phosphate buffer; stabilizer and preservative. (INCLUDED IN THE Waaler Rose KIT)

**QUALITY CONTROL**  
The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the controls below:

CONTROL (+) (INCLUDED IN THE Waaler Rose KIT)	<b>REF</b>	30.006.00
CONTROL (-) (INCLUDED IN THE Waaler Rose KIT)		30.006.00

**NECESSARY EQUIPMENT FOR TESTING**

- Reaction plate and disposable rod;
- Pipette, clock, or stopwatch.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**  
**A) TEST PROCEDURE**

**Qualitative Method**

1. Allow the R1, controls and samples to reach room temperature. Homogenize the R1 before its use.
2. Use a different circle of the slide for each sample or control.
3. Add 50 µL of the sample/control and 50 µL of the R1.
4. Homogenize it with a disposable rod and spread it all over the circle's surface. Use different rods for each sample.
5. Let the slide rest on a flat surface for 2 minutes. Then, gently tilt the slide at an 45° angle from the horizontal, returning to the starting position. Let it rest for 1 minute.
6. Immediately observe the result to avoid false positives. Avoid moving the slide during the observation. The visible presence of agglutination indicates a RF concentration equal or bigger than 8 UI/mL.

**Semi-quantitative method**

1. Dilute the sample with saline (0.9% NaCl) as follows:

Dilution	Saline (µL)	Sample(µL)	Dilution Factor
1/2	100	100 (pure)	2
1/4	100	100 (sol. 1/2)	4
1/8	100	100 (sol. 1/4)	8
1/16	100	100 (sol. 1/8)	16

2) Carry out the test of each dilution according to the qualitative method. The titer is defined as the highest dilution with a positive result.

**B) CALCULATIONS**  
The approximate concentration of RF in the sample is obtained by the following equation:

$$8 \times \text{Dilution Factor} = \text{UI/mL}$$

**C) INTERPRETATION**  
Rheumatoid factors (RF) are a group of antibodies that bind to the Fc fraction of immunoglobulin G. RF is frequently evaluated in patients with suspect of any form of arthritis, although positive results may be due to other causes and negative results do not exclude the illness. Nevertheless, in combination with signs and symptoms, it can play a role in the diagnosis and prognosis of diseases. High levels of RF occur in a variety of rheumatoid diseases, being present in 80% of rheumatoid arthritis cases and in 70% of Sjögren's syndrome cases. In addition, elevated RF levels are also found in: infections by the Epstein-Barr virus or Parvovirus infection; in 5-10% of healthy people (especially the elderly); chronic hepatitis; primary biliary cirrhosis; any chronic viral infection; bacterial endocarditis; leukemia; dermatomyositis; infectious mononucleosis; systemic sclerosis and systemic lupus erythematosus.

**PERFORMANCE CHARACTERISTICS**

Sensitivity
Detection Limit
8 IU/mL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
1000 mg/dL	20 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	84
Positive Concordance	100%
Negative Concordance	93,6%

**Prozone High Dose Effect:**  
There is no prozone effect up to 800 IU/mL of RF.

**Precision:**  
**Intra-assay:** the study was carried out in an analytical run with four replicates, for 10 days, using two positive and two negative samples, with 100% concordance.

**Inter-assay:** the study was carried out in an analytical run of four replicates, for 5 days, using two positive and two negative samples, with 100% concordance.

**RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

**REFERENCE RANGES**

Adults	Up to 8 UI/mL
--------	---------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

**WARNINGS AND PRECAUTIONS**

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechica.ind.br](http://www.biotechica.ind.br) or calling +55 35 3214 4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific

Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).

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- ASSIMEH, S. N.; JOHNSON, P. M. A Haemagglutination Method for Detection of Rheumatoid Factor Using Preserved Erythrocytes Covalently Coated With Human IgG. *Journal Of Immunological Methods*, v. 34, p.205-215, 1980.
- DORNER, Robert W.; ALEXANDER JUNIOR, Roy L; MOORE, Terry L. Rheumatoid Factors. *Clinica Chimica Acta*, v. 167, p.1-21, 1987.
- FULFORD, Karen M.; TAYLOR, Roger N.; PRZYBYSZEWski, Vncent A. Reference Preparation to Standardize Results of Serological Tests for Rheumatoid Factor. *Journal Of Clinical Microbiology*, v. 7, n° 5, p.434-441, 1978.
- KORITZ, T. N. et al. Measure of Rheumatoid Factor in Human Sera by Passive Haemagglutination of Human Erythrocytes Carrying Immunoglobulin Linked by Chromic Chloride. *Journal Of Immunological Methods*, v. 32, p.1-9, 1980.
- ROSE, Harry M. et al. Differential Agglutination of Normal and Sensitized Sheep Erythrocytes by Sera of Patients with Rheumatoid Arthritis. *Proc. Soc. Exp. Biol. Med.* v. 68, n. 1, p.1-6, 1948.
- SHMERLING, Robert H.; DELBANCO, Thomas L. The Reumathoid Factor: An Analysis of Clinical Utility. *The American Journal Of Medicine*, v. 91, p.528-534, 1991.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- WOLFE, Frederick; CATHEY, Mary Ann; ROBERTS, F. Kay. The Latex Test Revisited: Rheumatoid Factor Testing in 8,287 Rheumatic Disease Patients. *Arthritis And Rheumatism*, v. 34, n° 8, p.951-960, 1991.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: **AACC Press**, 2000

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