## Zinco

Zinc / Zinc  
Ref. 12.014.00 / FD 38.074.00 / SD 39.071.00

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#### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

#### INTENDED USE

Kit intended to determine zinc in serum and plasma samples. Diagnostic use only.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 30 days.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Colorimetric

Zinc ions (Zn<sup>2+</sup>) form with nitro-PAPS a colored complex that can be measured at 570-582 nm. The intensity of the formed color is proportional to the amount of zinc present in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum and plasma (heparin)

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Sample Stability
Serum and Plasma	2 to 8 °C	8 days
	-20 °C	1 year

#### PRODUCT DESCRIPTION

**R 1** Stabilizers; Surfactant and Preservative.

**R 2** Nitro-PAPS, Preservative.

**STD** Zinc Nitrate in solution. Traceable to reference material NIST-SRM 3168a.

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the controls below:

Pathological Control – Quantialt **REF** 13.004.00

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 578 nm (570 - 582 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

- Before starting, allow the reagents to reach a temperature of 15 to 30 ° C.
- Pipette into test tubes:

	Blank	STD	Sample
STD	-	50 µL	-
Sample	-	-	50 µL

	Purified water	50 µL	-	-
R1	800 µL	800 µL	800 µL	800 µL
Homogenize and incubate at 37 °C for 5 minutes				
R2	200 µL	200 µL	200 µL	200 µL

- Homogenize and incubate at 37 °C for 5 minutes
- Measure the absorbance of the sample and STD against the blank at 578 nm. The final color of the reaction is stable for 18 minutes.

#### B) CALCULATIONS

Zinc (µg/dL) =  $\frac{\text{Sample Absorbance}}{\text{Absorbance of STD}} \times \text{Concentration of STD (µg/dL)}$

**Exemplo:**

Concentration of STD = 200 µg/dL  
Sample Absorbance = 0,491  
Absorbance of STD = 0,263

Zinc (µg/dL) =  $\frac{0,491}{0,263} \times 200 = 373,4 \mu\text{g/dL}$

#### With Calibration Factor:

Calibration Factor =  $\frac{\text{Concentration of STD (µg/dL)}}{\text{Absorbance of STD}}$

Zinc (µg/dL) = Sample Absorbance x Calibration Factor

**Exemplo:**

Calibration Factor =  $\frac{200}{0,263} = 760,5$

Zinc (µg/dL) = 0,491 x 760,5 = 373,4 µg/dL

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicalda.ind.br](http://www.biotechnicalda.ind.br)

#### C) INTERPRETATION

Zinc is a cofactor in more than 70 enzymatic reactions, involving different metabolic processes such as the synthesis or degradation of carbohydrates, lipids, proteins and nucleic acids. Zinc deficiency can cause anemia, hepatosplenomegaly, slow growth and healing of ulcers and wounds, as well as changes in smell and taste. Low levels of zinc can occur in several physiological situations: pregnancy, contraceptive use and in some pathologies such as myocardial infarction, cirrhosis, lung infection, lymphoma and carcinoma.

#### PERFORMANCE CHARACTERISTICS

Operating range
47,54 to 819,73 µg/dL

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1,70 µg/dL	33,50 µg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Lipemia
25 mg/dL	20 mg/dL	Interfere

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 1.019x + 0.538
Correlation Coefficient (R)	0.9963

By applying the regression equation, the total systematic error estimated for decision levels of 80 µg/dL and 200 µg/dL was 2,57% and 2,17%, respectively.

**Precision:**

Determined with two runs in duplicate per day for 20 days.

Samples (µg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (µg/dL)	CV (%)	SD (µg/dL)	CV (%)
84,087	80	2,116	2,5	2,680	3,2
147,830	80	1,834	1,2	2,687	1,8
219,329	80	1,669	0,8	2,872	1,3

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity

has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.

- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES
80 - 120 µg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

- WARNINGS AND PRECAUTIONS**
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
  - The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling +55 35 3214 4646

- QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**
- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).
  - To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotecnica.ind.br](mailto:sac@biotecnica.ind.br).

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- Commission Directive 1999/11/EC of 8 March 1999 adapting to technical progress the principles of good laboratory practice as specified in Council Directive 87/18/EEC on the harmonisation of laws, regulations and administrative provisions relating to the application of the principles of good laboratory practice and the verification of their applications for tests on chemical substances (OJ L 77, 23.03.1999).
- MAKINO, T. A sensitive, direct colorimetric assay of serum zinc using Nitro-PAPS and microwell plates. *Clinica Chimica Acta.* n. 197, p. 209-220, 1991.
- GIROUX, E.I.; SCHECHTER, P. A Study of Zinc Distribution in Human Serum. *Bioinorganic Chemistry.* n. 5, p. 211-218, 1976
- WESTGARD, J. O. *et al.* A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v. 27, p.493-501, 1981.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests.* vol. 2, 5 ed. Washington DC: AACC Press, 2000.