Bilirrubina Total VET

Total Bilirubin VET / Bilirrubina Total VET
Ref. 90.053.00

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ENGLISH

INTENDED USE

Kit intended to determine total bilirubin in serum and plasma (heparin and EDTA) samples. Veterinary Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: DCA (Diazotized Dichloroaniline).

In acidic pH, direct and indirect bilirubin are coupled with diazotized dichloroaniline, being solubilized and producing a colored complex, called azobilirubin, that can be measured at 550 nm (505-580 nm). The color intensity is proportional to the concentration of total bilirubin in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (heparin and EDTA).

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	7 days
	-20 °C	6 months

*Necessary to keep out of the light when stored for longer than 8 hours.

PRODUCT DESCRIPTION

R 1	Surfactants, stabilizers, catalyst and detergent
R 2	Diazotized dichloroaniline < 2,0 mmol/L and stabilizer. pH < 5

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal VET	90.039.00
Normal Control – Qualitorm VET	90.040.00
Pathological Control – Quantialt VET	90.041.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 550 nm (505 – 580 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Sample Blank	Sample	Calibrator Blank	Calibrator
Calibrador	---	---	---	50 µL	50 µL
Sample	---	50 µL	50 µL	---	---
Purified water	50 µL	---	---	---	---
R1	1,0 mL	1,0 mL	1,0 mL	1,0 mL	1,0 mL
Homogenize and incubate the test tubes for 1 minute at 37°C					
R2	250 µL	---	250 µL	---	250 µL

- Homogenize and incubate at 37 °C for 4 minutes.
- Measure the Blank's, Sample Blank's, Sample's and Standard's absorbance against the purified water at 550 nm (505-580 nm). The final reaction is stable for 30 minutes.

B) CALCULATIONS

Total Bilirubin (mg/dL) = $\frac{\text{Sample's } \Delta A - \text{Blank}}{\text{Calibrator's } \Delta A - \text{Blank}}$ x Calibrator Concentration (mg/dL)

In which:

ΔA = Sample's or Calibrator's Absorbance – Sample or Calibrator Blank's Absorbance

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{Calibrator Concentration (mg/dL)}}{\text{Calibrator's } \Delta A - \text{Blank}}$

Total Bilirubin (mg/dL) = (Sample's ΔA – Blank) x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaldta.ind.br

C) INTERPRETATION

Assessing bilirubin levels is critical to understanding liver health, allowing for the diagnosis of conditions such as hepatitis and anemia. Bilirubin is processed in the liver and excreted in the bile, being the result of the breakdown of hemoglobin in the reticuloendothelial system. Its elevated levels may indicate various pathological conditions, including jaundice, which manifests itself when blood bilirubin is high, resulting in yellowish discoloration of the skin and mucous membranes due to the accumulation of bile pigments. Bilirubin measurement is useful in diagnosing diseases related to the liver, bile ducts, and hemolytic conditions. There are hereditary or acquired conditions that alter different stages of the process of production, uptake, storage, metabolism and elimination of bilirubin, affecting the balance between unconjugated (indirect), conjugated (direct) or total bilirubin, leading to hyperbilirubinemia. In newborns, the most common increase in indirect bilirubin is seen in physiological jaundice, which can trigger encephalopathy (kernicterus). Increased direct bilirubin is associated with inherited liver conditions such as Dubin-Johnson and Rotor diseases, damage to hepatocytes (caused by viruses, toxins, or alcohol), obstruction of bile ducts (stones or tumors), acute or chronic hepatitis, and adverse reactions to various drugs (such as chlorpromazine). On the other hand, the increase in indirect bilirubin is related to hemolytic anemias, autoimmune hemolysis, blood transfusions, resorption of hematomas, ineffective production of red blood cells, hereditary diseases such as Gilbert and Crigler-Najar diseases, in addition to neonatal jaundice.

PERFORMANCE CHARACTERISTICS

Operating range
0.12 to 50.00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.07 mg/dL	0.12 mg/dL

Analytical Specificity		
Species	Hemoglobin	Triglycerides
Canines	100 mg/dL	350 mg/dL
Felines	200 mg/dL	350 mg/dL
Equines	200 mg/dL	500 mg/dL
Bovines	200 mg/dL	350 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	40 in duplicate
Regression Equation	y = 0,943x + 0,0091
Correlation Coefficient (R)	0.9916
Accuracy - Felines	
Number of Samples	38 in duplicate
Regression Equation	y = 0.9448x - 0.2301
Correlation Coefficient (R)	0.9686
Accuracy - Equines	
Number of Samples	17 in duplicate
Regression Equation	y = 1.0296x - 0.0057
Correlation Coefficient (R)	0.9888
Accuracy - Bovines	
Number of Samples	32 in duplicate
Regression Equation	y = 0.994x + 0.0173
Correlation Coefficient (R)	0.9919

Precision:

Determined with two runs.

Samples (mg/dL)	Within-Run Precision			Samples (mg/dL)	Total Precision		
	SD (mg/dL)	CV (%)	CV (%)		SD (mg/dL)	CV (%)	CV (%)
0.31	0.01	4,7	0,31	0.01	0.01	4,6	
5.58	0.11	1,9	5.74	0.22	0.22	3,8	

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Canines	0,1 – 0,5 mg/dL
Felines	0,15 – 0,5 mg/dL
Equines	1,0 – 2,0 mg/dL
Bovines	0,01 – 0,5 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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